

**OneStep  
 Barbiturates Urine  
 RapiDip™ InstaTest**

REF 121062-1-19



For medical and other professional in vitro diagnostic use only.

**INTENDED USE**

OneStep Barbiturates Urine RapiDip™ InstaTest is a qualitative immunoassay intended to provide qualitative screening results for barbiturates in human urine at a cutoff concentration of 200ng/ml (Secobarbital). It is for health care professional use only.

*This assay provides only a preliminary result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography / Mass Spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are obtained.*

**SUMMARY AND EXPLANATION OF THE TEST**

Barbiturates are central nervous system depressants and used as hypnotic sedatives. Overdose and extended usage of barbiturates may lead to severe and/or permanent damage to the human nervous system. Barbiturates are classified as (1) ultra-short, (2) short-intermediate, and (3) long-acting. The duration range of the ultra short-acting compounds, secobarbital, pentobarbital etc. is from fifteen (15) minutes to six (6) hours. The duration range of the intermediate acting compounds, amobarbital, etc. is from three (3) to twenty-four (24) hours. The duration range of the long-acting compounds, phenobarbital etc. is from fifteen (15) to forty-eight (48) hours.

The most commonly abused barbiturates are short- and intermediate-acting agents. The long-acting agents are rarely subject to abuse. Barbiturate derivatives are excreted into urine in varying amounts of unchanged drug and metabolites. Long-acting barbiturates are excreted

with a higher percentage of unchanged drug in the urine, while shorter-acting barbiturates, secobarbital and amobarbital, are extensively metabolized and excreted in the urine with a smaller percentage of unchanged drugs.

OneStep Barbiturates Urine RapiDip™ InstaTest is designed to detect unchanged secobarbital in the urine; however, as with some other analytical methods such as EMIT and RIA, this assay can also detect other commonly encountered barbiturates, depending on the concentration of drug present in the sample. Phenobarbital positives have been noted in chronic users up to several weeks after cessation of use. With standard single doses of secobarbital, pentobarbital, or amobarbital, positive results may be identified from 30 hours to 76 hours.

**TEST PRINCIPLE**

This assay is a one-step lateral flow chromatographic immunoassay. The test strip includes 1) a burgundy-colored conjugate pad containing mouse anti-barbiturate antibodies coupled to colloidal gold; and 2) nitrocellulose membrane containing a Test (T) line and a Control (C) line. The Test line is coated with barbiturate-BTG, and the Control line is coated with goat anti-rabbit IgG antibody.

OneStep Barbiturates Urine RapiDip™ InstaTest is a competitive binding immunoassay. The barbiturate in the urine specimen competes with the barbiturate-BTG antigen coated on the nitrocellulose membrane for the limited binding sites of the conjugated anti-barbiturate antibodies.

When an adequate amount of urine specimen is applied to the sample pad of the device, the urine specimen migrates by capillary action through the test strip. If the level of barbiturate in the urine specimen is below the cutoff (200 ng/ml), the Test line should appear as a visible burgundy line. If the level of barbiturate in the urine specimen is at or above the cutoff, no Test line develops.

The C line binds to the gold-conjugated rabbit IgG and forms a burgundy color line, regardless of the presence of barbiturate.

**SPECIMEN COLLECTION AND PREPARATION**

1. Each urine specimen must be collected in a clean container. Do not mix specimens.
2. Specimens may be kept at 15-30°C (59-86°F) for 8 hours, at 2-8°C for up to 3 days and at -20°C or lower for long term storage.

**MATERIALS AND COMPONENTS**

Diagnostic Automation/ Cortez Diagnostics, Inc.  
 21250 Califa St, Suite 102 and 116, Woodland Hills, CA 91367 USA Phone: 818-591-3030, Fax : 818-591-8383  
 Email: onestep@rapidtest.com Website: www.rapidtest.com

**Materials provided with the test kits**

- 50 test strips each sealed in a pouch with desiccant.
- 1 package insert (Instructions for Use).

**Materials required but not provided**


- Specimen collection containers
- Timer

**ASSAY PROCEDURE**

1. Refrigerated specimens and other test materials, including devices, **must be equilibrated to room temperature before testing.**
2. Open the foil pouch at the notch and remove the test device.
3. Dip the test device in the specimen for at least 10 seconds. Keep the specimen surface at the level indicated by the arrow sign on the device.
4. Remove the test device from the specimen, and place it on a flat, dry surface.
5. Read the test result between four (4) to seven (7) minutes after adding the specimen.

**RESULTS**

**IMPORTANT: Do not read test results after seven (7) minutes. The T Line should always be interpreted independently of the C Line.**



**Positive:**

If only the C line appears, the test indicates that the level of barbiturates in the specimen is at or above the cutoff (200ng/ml).

**Samples with preliminary positive results should be confirmed with a more specific method before a positive conclusion is made.**

**Negative:**

If both C line and T line appear, the test indicates that the level of barbiturates in the specimen is below the cutoff (200ng/ml).

**Note: A very faint T line should be considered negative.**

**Invalid:**

If no C line develops within 5 minutes, repeat the assay with a new test device.

**PERFORMANCE CHARACTERISTICS**

**1. Accuracy**

A study was performed at three different Physician's Office Laboratories (POL) and a Reference Laboratory. One hundred (100) clinical samples were blind labeled and tested. Each sample was tested at each site, and compared with GC/MS results.

The results agreed 100% with the GC/MS data of specimens at levels below the cutoff (negative) and above 125% of the cutoff (positive). One (1) discrepancies was observed on the specimens at the level between the cutoff and 125% of the cutoff.

The overall agreement was 99.8%.

		BAR Test		Total	Agreement
		Positive	Negative		
GC/MS (ng/ml)	Drug-free	0	200	200	100%
	<75% (0-150)	0	12	12	100%
	75%-Cutoff (150-200)	0	20	20	100%
	Cutoff-125% (200-250)	27	1	28	96.4%
	Positive (>250)	140	0	140	100%
<b>Total</b>		167	233	400	99.8%

**2. Precision**

The precision was determined by replicate assays of four different levels of samples with three different production lots. The devices were tested for five consecutive days five times each, for a total of 25 assays for each control.

The results indicate 100% precision for the replicate within each lot and no appreciable inter-lot variation occurred across the three different lots of devices.

**3. Cross-Reactivity**

To determine the cross-reactivity of the structurally related compounds with the device, the following compounds were spiked into known drug-free urine pools and tested. Those compounds showed a positive response at the concentration indicated in the following table:

Description	Concentration (ng/ml)
Amobarbital	250
Barbital	250
Butobarbital	300
Butalbital	200
Phenobarbital	200
Pentobarbital	250
Secobarbital	200

**4. Interference**

To determine the interference of structurally unrelated analytes, the following analytes were spiked into known drug-free urine pools, as well as the Secobarbital positive (spiked with Secobarbital to the level of 200 ng/ml) urine pools and were tested. No significant interference with either negative or positive results was observed at the concentrations listed in the following table:

**Compounds listed in this table found not to interfere with the test results at the concentration of 1 mg/ml:**

Acetaminophen	Cortisone
Acetylsalicylic Acid	Dextromethorphan
Amikacin	Ethanol
Amitriptyline	Lidocaine
Ampicillin	Methadone
Arterenal	Methanol
Aspirin	Oxalic Acid
Atropine	Penicillin-G (Benzylpenicillin)
Benzoic Acid	Pheniramine
Benzoylecgonine	Phenylpropanolamine
Caffeine	Ranitidine
(+)-Chlorpheniramine	Salicylic Acid
Cocaine	Thioridazine
Codeine	Trifluoperazine

Biological Analytes	Concentration
Albumin	2 mg/ml
Bilirubin	1 mg/ml
Creatine	1 mg/ml
Hemoglobin	1 mg/ml
Glucose	2 mg/ml
pH	5.0 – 9.0
Vitamin C (L-Ascorbic Acid)	1 mg/ml
Uric Acid	1 mg/ml

*There is a possibility that other substances and/or factors not listed, may interfere with the test and cause false results.*

**STORAGE**

Store the kit at room temperature 15-30°C (59-86°F). Each device may be used until the expiration date printed on the label if it remains sealed in its foil pouch containing desiccant.

**Do not freeze and/or expose the kit to temperatures over 30 °C (86°F).** 15°C  30°C

**QUALITY CONTROL**

- Built-in Control Features**  
This test contains a built-in control feature, the C line. The presence of the C line indicates that an adequate sample volume was used and that the reagents migrated properly. If a C line does not form, the test is considered invalid. In this case, review the whole procedure and repeat the testing with a new device.
- External Quality Control**  
Users should always follow the appropriate federal, state, and local guidelines concerning the running of external quality controls. SAMHSA<sup>3</sup> recommends that the concentration of drug(s) in positive and negative controls be approximately 25% above and below the cutoff concentration of the assay.

**LIMITATIONS OF PROCEDURE**

- OneStep Barbiturates Urine RapiDip™ InstaTest is for professional *in vitro* diagnostic use only.
- Results obtained by this device provide only a preliminary qualitative result. A more specific alternate chemical method must be used in order to obtain a confirmed result.
- This product is designed for testing human urine only.
- Adulterants such as bleach or other strong oxidizing agents may produce erroneous test results if present in the sample. When

suspected, collect a fresh specimen and repeat the test with a new device.

5. Samples in which bacterial contamination is suspected should not be used. These samples may interfere with the test and cause false results.

**EXPECTED VALUES**

OneStep Barbiturates Urine RapiDip™ InstaTest is designed to detect barbiturates (secobarbital) in human urine at a cutoff concentration of 200ng/ml.

**PRECAUTION**

1. The instructions must be followed exactly to obtain accurate results.
2. Do not open the sealed pouch, unless ready to conduct the assay.
3. Do not use expired devices.
4. Dispose of all specimens and used assay materials as potentially biohazardous.

**REFERENCE**

1. FDA Guidance for Labeling Urine Drugs of Abuse Screening Testing, Kshit Mohan, 7/21/87.
2. Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, Fed. Register. 53 (69): 11970 (1988).
3. Urine Testing for Drugs of Abuse, National Institute on Drug Abuse (NIDA): Research Monograph 73, 1986.
4. Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man, 4th ED., Biomedical Publ., Davis, CA; pp 68-69, 1995.
5. Wilson, John, Abused Drugs II, a Laboratory Pocket Guide., AACC Press., Washington, DC; 1994.

<p><b>ISO 13485</b> <b>ISO 9001</b></p> 	
 <p><b>Diagnostic Automation/ Cortez Diagnostics, Inc.</b>  <b>21250 Califa St, Suite 102 and 116, Woodland Hills, California 91367 USA</b></p>	
<b>Date Adopted</b>	<b>2016-10-28</b>
REF 121062-1-19	<p><b>OneStep Barbiturates Urine RapiDip™ InstaTest</b></p>
EC REP	<p><b>CEpartner4U, Esdoornlaan 13, 3951DB Maarn. The Netherlands.</b>  <a href="http://www.cepartner4u.eu">www.cepartner4u.eu</a></p>
<p>Revision D Evision Date: 2007-01-29</p>	

Diagnostic Automation/ Cortez Diagnostics, Inc.  
 21250 Califa St, Suite 102 and 116, Woodland Hills, CA 91367 USA Phone: 818-591-3030, Fax : 818-591-8383  
 Email: onestep@rapidtest.com Website: www.rapidtest.com