

**AccuDiag™
EBV-EA IgG
ELISA Kit**

REF 1415-2



Test	EBV-EA IgG ELISA
Method	Enzyme Linked Immunosorbent Assay
Principle	Qualitative
Sample	10 µL serum
Total Time	~ 65 min.
Shelf Life	12 Months from the manufacturing date
Sensitivity	81 %
Specificity	92 %

INTENDED USE

The Diagnostic Automation, Inc. EBV-EA IgG ELISA test system provides a means for the qualitative detection of IgG antibodies to Epstein-Barr Virus (EBV) Early Antigen (EA) in human sera. This test system is intended to aid in the diagnosis of infectious mononucleosis (IM) when used with other EBV serological tests. The performance characteristics have not been established to aid in the diagnosis of acute IM. The test is intended for *in vitro* diagnostic use only.

SUMMARY AND EXPLANATION

Epstein - Barr virus (EBV) causes infectious mononucleosis; a self-limiting lymphoproliferative disease. EBV is a ubiquitous human virus. By adulthood virtually everyone has been infected and has developed immunity to the virus. In underdeveloped countries, seroconversion to the virus takes place in early childhood and is usually asymptomatic. In more affluent countries, primary EBV infections are often delayed until adolescence or later, and manifest as IM in about 50% of this age group.

Following seroconversion, whether symptomatic or not, EBV establishes a chronic latent infection in B lymphocytes which lasts probably for life. EBV replicates in oropharyngeal epithelial cells and is present in the saliva of most patients with IM. Also, 10-20% of healthy persons who are EBV antibody positive shed the virus in their oral secretions. Reactivation of the latent viral carrier state, as evidenced by increased rates of virus shedding, is enhanced by immunosuppression, pregnancy, malnutrition or disease. Chronic EBV infections, whether latent or active, are rarely associated with disease. However, EBV has been implicated at least a contributing factor in the etiology of nasopharyngeal carcinoma, Burkitt's lymphoma, and lymphomas in immunodeficient patients.

The Paul-Bunnell-Davidsohn test for heterophile antibody is highly specific for IM. However, 10-15% of adults, and higher percentage of children and infants with primary EBV infections do not develop heterophile antibodies. EBV

specific tests are needed to differentiate primary EBV infections that are heterophile negative from mononucleosis-like illnesses caused by other agents, such as cytomegalovirus, adenovirus, and Toxoplasma gondii.

Antibody titers to specific EBV antigens correlate with different stages of IM. Both IgM and IgG antibodies to the viral capsid antigen (VCA) peak 3 to 4 weeks after primary EBV infection. IgM anti VCA declines rapidly and is usually undetectable after 12 weeks. IgG anti-VCA titers decline slowly after peaking but last indefinitely. Antibodies to EBV nuclear antigen (EBV-NA) detected by anticomplement immunofluorescence develop from 1 month to 6 months after infection; and, like anti-VCA, persist indefinitely. Antibodies to EBV-NA indicate that the EBV infection was not recent.

EBV early antigen (EA) consists of two components; diffuse (D), and restricted (R). The terms D and R reflect the different patterns of immunofluorescence staining exhibited by the two components. Antibodies to EA may appear transiently for up to three months or longer during the acute phase of IM in 85 % of patients. The antibody response to EA in IM patients is usually to the D component, whereas silent seroconversion to EBV in children may produce antibodies to the R components. A definitive diagnosis of primary EBV Infection can be made with 95 % of acute phase sera based on antibody titers to VCA, EBV-NA, and EA.

Antibodies to EA, usually to the R component, together with antibodies to EBV-NA and high titers of IgG anti-VCA, may be associated with reactivation of the latent viral carrier state. EBV positive serology associated with reactivation of EBV is found in sera of patients with immunodeficiencies, patients with recurrent parotitis, immunosuppressed patients, pregnant women, and persons of advanced age. Antibodies to the R component may be found at moderate to high levels in patients with Burkitt's lymphoma. In contrast patients with nasopharyngeal carcinoma may produce high titer antibodies to the D component.

Elevated levels of anti-EA and IgG anti-VCA may be detected in patients with chronic or recurrent illness suspected of being caused by EBV. However, a diagnosis of chronic EBV should not be based on the presence of antibodies to EA since elevated anti-EA titers may also be found in patients with other diseases as well as in healthy individuals with past EBV infections.

TEST PRINCIPLE

The Diagnostic Automation, Inc. EBV-EA ELISA Test system is designed to detect IgG class antibodies to Epstein - Barr virus Early Antigen in human sera. Wells of plastic microwell strips are sensitized by passive absorption with EBV-EA antigen. The test procedure involves three incubation steps:

1. Test sera (properly diluted) are incubated in antigen coated microwells. Any antigen specific antibody in the sample will bind to the immobilized antigen .The plate is washed to remove unbound antibody and other serum components.
2. Peroxidase Conjugated goat anti-human IgG is added to the wells and plate is incubated. The Conjugate will react with IgG antibody immobilized on the solid phase in step 1. The wells are washed to remove unreacted Conjugate.
3. The microwells containing immobilized peroxidase Conjugate are incubated with peroxidase Substrate Solution. Hydrolysis of the Substrate by peroxidase produces a color change. After a period of time the reaction is stopped and the color intensity of the solution is measured photometrically. The color intensity of the solution depends upon the antibody concentration in the original test sample.

SPECIMEN COLLECTION AND PREPARATION

1. It is recommended that specimen collection be carried out in accordance with CLSI document M29: Protection of Laboratory Workers from infectious Disease (Current Edition).



2. No known test method can offer complete assurance that human blood samples will not transmit infection. Therefore, all blood derivatives should be considered potentially infectious.
3. Only freshly drawn and properly refrigerated sera obtained by approved aseptic venipuncture procedures should be used in this assay. No anticoagulants or preservatives should be added. Avoid using hemolyzed, lipemic, or bacterially contaminated sera.
4. Store sample at room temperature for no longer than 8 hours. If testing is not performed within 8 hours, sera may be stored between 2 °C and 8°C for no longer than 48 hours. If delay in testing is anticipated, store test sera at -20°C or lower. Avoid multiple freeze/thaw cycles that may cause loss of antibody activity and give erroneous results. It is the responsibility of the individual laboratory to use all available references and/or its own studies to determine stability criteria for its laboratory (32).

12. Disposal basin and disinfectant (example: 10 % household bleach - 0.5 % sodium hypochlorite).

ASSAY PROCEDURE

1. Remove the individual components from storage and allow them to warm to room temperature (20-25°C).
2. Determine the number of microwells needed. Allow six Control/Calibrator determinations (one Blank, one Negative Control, three Calibrators and one Positive Control) per run. A reagent Blank should be run on each assay. Check software and reader requirements for the correct Controls/Calibrator configurations. Return unused strips to the resealable pouch with desiccant, seal, and return to storage between 2 ° and 8 °C.

MATERIALS AND COMPONENTS

Materials provided with the test kits

Each Test System contains the following components in sufficient quantities to perform the number of tests indicated on packaging label. **NOTE: The following components contain Sodium Azide as a preservative at a concentration of < 0.1 % (w/v): Controls, Calibrators and Sample Diluent.**

1. **Plate:** 96 wells configured in twelve 1x8 well strips coated with inactivated EBV-EA antigen derived from induced Raji cells. The strips are packaged in a strip holder and sealed in an envelope with desiccant.
2. **Conjugate:** Conjugated (horseradish peroxidase) goat anti-human IgG (Fc chain specific). One, 15 mL vial with a white-capped bottle. Ready to use.
3. **Positive Control (Human Serum):** One, 0.35mL vial with a red capped.
4. **Calibrator (Human Serum):** One 0.5 mL vial with a blue cap.
5. **Negative Control (Human Serum):** One, 0.35mL vial with a green capped
6. **Sample Diluent:** One, 30mL, green-capped, bottle containing Tween-20, bovine serum albumin and phosphate-buffered-saline (pH 7.2 ± 0.2). Ready to use. **NOTE:** The sample diluents will change color when combined with serum.
7. **TMB:** One 15 mL amber bottle (amber cap) containing 3,3',5,5'-tetramethylbenzidine (TMB). Ready to use.
8. **Stop Solution:** One 15 mL bottle (red cap) containing 1 M H₂SO₄, 0.7 M HCl. Ready to use.
9. **Wash buffer Concentrate (10X):** Dilute 1 part concentrate + 9 parts deionized or distilled water. One 100 mL bottle (clear cap) containing a 10X concentrated phosphate-buffered-saline and Tween-20 solution (blue solution). **NOTE: 1 X solution will have a pH of 7.2 ± 0.2.**

The following components are not kit lot number dependent and may be used interchangeably with DAI ELISA assays: TMB, Stop Solution, and Wash Buffer.

Note: Kit also contains:

1. Component list containing lot specific information is inside the kit box.
2. Package inserts providing instructions for use.

Materials required but not provided

1. ELISA microwell reader capable of reading at a wavelength of 450 nm.
2. Pipettes capable of accurately delivering 10 to 200 µL.
3. Multichannel pipette capable of accurately delivering 50-200 µL.
4. Reagent reservoirs for multichannel pipettes.
5. Wash bottle or microwell washing system.
6. Distilled or deionized water.
7. One liter graduated cylinder.
8. Serological pipettes.
9. Disposable pipette tips.
10. Paper towels.
11. Laboratory timer to monitor incubation steps.

EXAMPLE PLATE SET-UP		
	1	2
A	Blank	Patient 3
B	Neg. Control	Patient 4
C	Calibrator	Etc.
D	Calibrator	
E	Calibrator	
F	Pos. Control	
G	Patient 1	
H	Patient 2	

3. Prepare a 1:21 dilution (e.g.: 10 µL of serum + 200 µL of Sample Diluent.) of the Negative Control, Calibrator, Positive Control, and each patient serum. **NOTE: The Sample Diluent will undergo a color change confirming that the specimen has been combined with the diluent.**
4. To individual wells, add 100 µL of each diluted control, calibrator and sample. Ensure that the samples are properly mixed. Use a different pipette tip for each sample.
5. Add 100 µL of Sample Diluent to well A1 as a reagent blank. Check software and reader requirements for the correct reagent blank well configuration.
6. Incubate the plate at room temperature (20-25°C) for 25 ± 5 minutes.
7. Wash the microwell strips 5 X.
 - A. **Wash the microwell Procedure:**
 - a. Vigorously shake out the liquid from the wells.
 - b. Fill each microwell with Wash Buffer. Make sure no air bubbles are trapped in the wells.
 - c. Repeat steps a. and b. for a total of 5 washes.
 - d. Shake out the wash solution from all the wells. Invert the plate over a paper towel and tap firmly to remove any residual wash solution from the wells. Visually inspect the plate to ensure that no residual wash solution remains. Collect wash solution in a disposable basin and treat with disinfectant at the end of the days run.
 - B. **Automated Wash Procedure:**

If using automated microwell wash system, set the dispensing volume to 300-350 µL/well. Set the wash cycle for 5 washes with no delay between washes. If necessary, the microwell plate may be removed from the washer, inverted over a paper towel and tapped firmly to remove any residual wash solution from the microwells.
8. Add 100 µL of the Conjugate to each well, including reagent blank well, at the same rate and in the same order as the specimens were added.
9. Incubate the plate at room temperature (20-25°C) for 25 ± 5 minutes.
10. Wash the microwells by following the procedure as described in step 7.
11. Add 100 µL of TMB to each well, including reagent blank well, at the same rate and in the same order as the specimens were added.
12. Incubate the plate at room temperature (20-25°C) for 10 to 15 minutes.
13. Stop the reaction by adding 50 µL of Stop Solution to each well, including reagent blank well, at the same rate and in the same order as the TMB was added. Positive samples will turn from blue to yellow. After adding the Stop Solution, tap the plate several times to ensure that the samples are thoroughly mixed.



- Set the microwell reader to read at a wavelength of 450 nm and measure the optical density (OD) of each well against the reagent blank. The plate should be read within 30 minutes after the addition of the Stop Solution.

Calibrator	≥ 0.300
Positive Control	≥ 0.500

RESULTS

A. Calculations:

1. Correction Factor

A cutoff OD value for positive samples has been determined by the manufacturer and correlated to the Calibrator. The correction factor (CF) will allow you to determine the cutoff value for positive samples and to correct for slight day to day variations in test results. The correction factor is determined for each lot of kit components and is printed on the data label located in the kit box.

2. Cutoff OD Value

To obtain the cutoff OD value, multiply the CF by the mean OD of the Calibrator determined above.
 (CF x mean OD of Calibrator = cutoff OD value)

3. Index Values or OD Ratios

Calculate the Index Value or OD Ratio for each specimen by dividing its OD value by the cutoff OD from step 2.

Example:

Mean OD of Calibrator = 0.793
 Correction Factor (CF) = 0.25
 Cut off OD = 0.793 x 0.25 = 0.198
 Unknown Specimen OD = 0.432
 Specimen Index Value or OD Ratio = 0.432 / 0.198 = 2.18

B. Interpretations:

Index Values or OD ratios are interpreted as follows:

	Index Value or OD Ratio
Negative Specimens	≤ 0.90
Equivocal Specimens	0.91 to 1.09
Positive Specimens	≥ 1.10

- An OD ratio ≤ 0.90 indicates no detectable IgG antibodies to EBV-EA.
- An OD ratio ≥ 1.10 indicates that IgG antibodies specific to EBV-EA were detected.
- Specimens with OD ratio values in the equivocal range (0.91 - 1.09) should be retested in duplicate. Report any two of the three results which agree. Evaluate repeatedly equivocal specimen using an alternate serological method and/or re-evaluate by drawing another sample one to three weeks later.

QUALITY CONTROL

- Each time the assay is run the Calibrator must be run in triplicate. A reagent blank, Negative Control, and Positive Control must also be included in each assay.
- Calculate the mean of the three Calibrator wells. If any of the values differ by more than 15 % from the mean, discard that value and calculate the mean using the remaining two wells.
- The mean OD value for the Calibrator and the OD values for the Positive and Negative Controls should fall within the following ranges:

	OD Range
Negative Control	≤ 0.250

- The OD of the Negative Control divided by the mean OD of the Calibrator should be ≤ 0.9.
 - The OD of the Positive Control divided by the mean OD of the Calibrator should be ≥ 1.25.
 - If the above conditions are not met the test should be considered invalid and should be repeated.
- The Positive Control and Negative Control are intended to monitor for substantial reagent failure and will not ensure precision at the assay cut-off.
 - Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations.
 - Refer to CLSI document C24: Statistical Quality Control for Quantitative Measurements for guidance on appropriate QC practices.

PERFORMANCE CHARACTERISTICS

1. Comparative Study:

A comparative study was performed to demonstrate the equivalence of the Diagnostic Automation, Inc. EBV-EA IgG ELISA test system to the Diagnostic Automation, Inc. IFA EBV-EA test system when evaluating specimens submitted for the diagnosis of IM.

The performance of the Diagnostic Automation, Inc. EBV-EA ELISA test system was evaluated in a three-site clinical investigation. Briefly, there were a total of 273 specimens tested: 125 at site one, 125 at site two, and 23 at site three. Clinical specimens tested at sites one and two consisted of a mixture of routine specimens which were sent to a reference laboratory in Southwestern United States for normal EBV serological analysis, and normal donor specimens. Repository specimens were tested at site 3. These had been previously tested and were found to be positive for antibody to EBV-EA. Equivocal specimens were excluded from any further analysis. Table 1 below shows the results from Clinical Site # 1. Table 2 shows the results from Clinically Site # 2. Results from Clinical Site # 3 have not been shown separately, but are included in Table 3 which shows the results of all 3 clinical sites combined.

Table 1: Calculation of Relative Sensitivity, Specificity and Agreement; Study Site One.

DAI IFA EBV-EA IgG	DAI EBV-EA IgG Result	DAI EBV-EA IgG Result			Totals
		+	-	±	
+	+	22	9	2	33
-	-	2	73	3	78
±	±	5	5	4	14
Totals		29	87	9	125

Relative Sensitivity = 22/31 = 71% (95% Confidence Interval * = 55 to 87%)
 Relative Specificity = 73/75 = 97% (95 % Confidence Interval * = 94 to 100 %)
 Relative Agreement = 95/106 = 90 % (95 % Confidence Interval * = 84 to 95 %)

*95 % confidence Intervals calculated using the exact method.

Table 2: Calculations of Relative Sensitivity, Specificity, and Agreement; Study Site Two.

DAI IFA EBV-EA IgG	DAI EBV-EA IgG Result	DAI EBV-EA IgG Result			Totals
		+	-	±	
+	+	25	6	4	35
-	-	7	83	0	90
±	±	0	0	0	0
Totals		32	89	4	125

Relative Sensitivity = 25/31 = 81% (95% Confidence Interval * = 67 to 94%)
 Relative Specificity = 83/90 = 92% (95 % Confidence Interval * = 87 to 98 %)
 Relative Agreement = 108/121 = 89 % (95 % Confidence Interval* = 84 to 95 %)
 *95 % confidence Intervals calculated using the exact method.

Table 3: Calculation of Relative Sensitivity, Specificity and Agreement; All three sites combined.

DAI IFA	DAI EBV-EA IgG Result	DAI EBV-EA IgG Result			Totals
		+	-	±	
+	+	67	15	7	89
EBV-EA	-	11	156	3	170
	±	5	5	4	14
	Totals	83	176	14	273

Relative Sensitivity = 67/82 = 82% (95% Confidence Interval * = 67 to 94%)
 Relative Specificity = 156/167 = 93% (95 % Confidence Interval * = 90 to 97 %)
 Relative Agreement = 223/249 = 90 % (95 % Confidence Interval* = 86 to 93 %)
 *95 % confidence Intervals calculated using the exact method.

In table 3 above, there were a total of 26 discrepant samples. There were 11/273 (4.0%) samples which were positive on the ELISA and negative by IFA. There were a total of 15/273 (5.6%) samples which were negative on the ELISA and positive by IFA. Of this group of 15 samples 1/15 had an endpoint titer on IFA of 1:160, 4/15 samples had an endpoint titer of 1:20, and the remainder (10/15) had an endpoint on IFA of 1:10.

NOTE: Be advised that relative refers to the comparison of this assay's results to that of a similar assay. There was not an attempt to correlate the assay's results with disease presence or absence. No judgment can be made on the comparison assay's accuracy to predict disease.

All 273 specimens tested during the clinical study were also tested for EBV-VCA IgG, EBV-VCA IgM, EBNA-1 IgG, and heterophile antibody. The results of this ancillary testing is depicted in Table 4.

Table 4. Results of the EBV-EA IgG ELISA by category of EBV serological status.

EBV-EA IgG ELISA Results		Acute	Seropositive	Seronegative
		VCA IgM Positive Heterophile Positive EBNA Negative	VCA IgG Positive EBNA Positive VCA IgM Negative Heterophile Negative	VCA IgG Negative VCA IgM Negative EBNA Negative Heterophile Negative
Positive		0	54	1
Negative		2	106	29
Equivocal		0	10	1
TOTALS		2/273	170/273	31/273

2. Reproducibility:

Reproducibility was evaluated as outlined in the FDA guidance document: review Criteria for In Vitro Diagnostic Devices for Detection of IgM antibodies to Viral Antigens. Reproducibility studies were conducted at both clinical sites using the same specimens.

Briefly, six specimens were tested, two relatively strong positive specimens, two specimens near the cut-off, and two which clearly negative. Additionally, the kit's negative control and positive control were included as additional panel members at site one, for a total of eight specimens. On each day of testing, each of the eight specimens were assayed in eight replicate wells. Testing was performed for a total of three days at each site. A summary of this investigation appears in Table 5 below:

Table 5: Summary of Precision Testing Conducted at Clinical Sites 1 and 2

Specimen	Site	Mean Ratio	Result	Standard Deviation	Overall % CV
1	1	3.46	Pos.	0.55	16.0
	2	3.06		0.31	10.0
2	1	2.78	Pos.	0.60	21.7
	2	2.76		0.30	10.9
3	1	1.81	Near Cutoff	0.34	18.7
	2	1.69		0.28	16.6
4	1	1.15	Near Cutoff	0.19	16.5
	2	1.13		0.26	23.0
5	1	0.27	Neg.	0.06	24.0
	2	0.37		0.07	17.8
6	1	0.28	Neg.	0.08	30.1
	2	0.75		0.10	13.5
Positive Control	1	3.71	Pos.	1.00	27.0
	2	3.49		0.47	13.4
Negative Control	1	0.28	Neg.	0.06	23.6
	2	0.28		0.03	12.2

3. Cross Reactivity:

A Study was performed to investigate the possibility of cross-reactivity with other viruses. In this study, ten specimens were evaluated. Two of the specimens were strongly reactive to Rubella, two to Rubeola, two to HSV-1, two to HSV-2, and two to CMV. None of the ten specimens were reactive on the EBV-EA IgG ELISA test system, indicating that there is little potential for cross-reactivity with such patient specimens.

LIMITATIONS OF PROCEDURE

1. A diagnosis should not be made on the basis of anti-EBV-EA results alone. Test results alone. Test results for anti-EBV-EA should be interpreted in conjunction with the clinical evaluation and the results of other diagnostic procedures.
2. This test detects both R and D components of EA. The test system is not designed to differentiate between antibodies to the R and D components.
3. The use of hemolytic, lipemic, bacterially contaminated or heat inactivated specimens should be avoided. Erroneous results may occur.
4. The assay performance characteristics have not been established for matrices other than sera.
5. The magnitude of the measured result above the cut-off is not indicative of the total amount of antibody present and cannot be correlated to IFA titers.
6. Assay Performance characteristics have not been established for visual result determinations.
7. Caution should be used when evaluating samples obtained from immunosuppressed patients.
8. Prevalence of the analyte will affect the assay's predictive value.
9. The performance characteristics of this assay have not been established for Burkitt's lymphoma, nasopharyngeal carcinoma, and lymphoproliferative disorders. The performance has been established for the aid in the diagnosis of EBV-associated infectious mononucleosis.



EXPECTED VALUES

To establish or estimate the expected reactivity rate, the 250 specimens which were tested at the two clinical sites were analyzed. This represented two groups of specimens; 150 clinical specimens which were sent to the laboratory for routine EBV serological analysis, and 100 random normal donor specimens. With respect to the clinical population, 43/150 (28.7%) were positive, 101/150 (67.3%) were negative, and 6/150 (4.0%) were equivocal. With respect to the normal population, 18/100 (18%) were positive, 76/100 (76%) were negative, and 6/100 (6%) were equivocal.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. Normal precautions exercised in handling laboratory reagents should be followed. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves, and eye/face protection. Do not breathe vapor. Dispose of waste observing all local, state, and federal laws.
3. The wells of the ELISA plate do not contain viable organisms. However, the strips should be considered POTENTIALLY BIOHAZARDOUS MATERIALS and handled accordingly.
4. The human serum controls are POTENTIALLY BIOHAZARDOUS MATERIALS. Source materials from which these products were derived were found negative for HIV-1 antigen, HBsAg, and for antibodies against HCV and HIV by approved test methods. However, since no test method can offer complete assurance that infectious agents are absent, these products should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institute of Health manual "Biosafety in Microbiological and Biomedical Laboratories": current edition; and OSHA's Standard for Bloodborne Pathogens.
5. Adherence to the specified time and temperature of incubations is essential for accurate results. All reagents must be allowed to reach room temperature (20-25°C) before starting the assay. Return unused reagents to refrigerated temperature immediately after use.
6. Improper washing could cause false positive or false negative results. Be sure to minimize the amount of any residual wash solution; (e.g., by blotting or aspiration) before adding Conjugate or Substrate. Do not allow the wells to dry out between incubations.
7. The Sample diluent, controls, wash buffer, and conjugate contain sodium azide at a concentration of 0.1% (w/v). Sodium azide has been reported to form lead or copper azides in laboratory plumbing which may cause explosions on hammering. To prevent, rinse sink thoroughly with water after disposing of solution containing sodium azide.
8. The Stop Solution is TOXIC. Causes burns. Toxic by inhalation, in contact with skin and if swallowed. In case of accident or if you feel unwell, seek medical advice immediately.
9. The TMB Solution is HARMFUL. Irritating to eyes, respiratory system and skin.
10. The Wash Buffer concentrate is an IRRITANT. Irritating to eyes, respiratory system and skin.
11. Wipe bottom of plate free of residual liquid and/or fingerprints that can alter optical density (OD) readings.
12. Dilution or adulteration of these reagents may generate erroneous results.
13. Reagents from other sources or manufacturers should not be used.
14. TMB Solution should be colorless, very pale yellow, very pale green, or very pale blue when used. Contamination of the TMB with conjugate or other oxidants will cause the solution to change color prematurely. Do not use the TMB if it is noticeably blue in color.
15. Never pipette with mouth. Avoid contact of reagents and patient specimens with skin and mucous membranes.
16. Avoid microbial contamination of reagents. Incorrect results may occur.

17. Cross contamination of reagents and/or samples could cause erroneous results.
18. Reusable glassware must be washed and thoroughly rinsed free of all detergents.
19. Avoid splashing or generation of aerosols.
20. Do not expose reagents to strong light during storage or incubation.
21. Allowing the microwell strips and holder to equilibrate to room temperature prior to opening the protective envelope will protect the wells from condensation.
22. Wash solution should be collected in a disposal basin. Treat the waste solution with 10% household bleach (0.5% sodium hypochlorite). Avoid exposure of reagents to bleach fumes.
23. Caution: Liquid waste at acid pH should be neutralized before adding to bleach solution.
24. Do not use ELISA plate if the indicator strip on the desiccant pouch has turned from blue to pink.
25. Do not allow the conjugate to come in contact with containers or instruments that may have previously contained a solution utilizing sodium azide as a preservative. Residual amounts of sodium azide may destroy the conjugate's enzymatic activity.
26. Do not expose any of the reactive reagents to bleach-containing solutions or to any strong odors from bleach-containing solutions. Trace amounts of bleach (sodium hypochlorite) may destroy the biological activity of many of the reactive reagents within this kit.

STORAGE

1. Store the unopened kit between 2 °C and 8 °C.
2. Coated microwell strips: Store between 2 ° and 8 °C. Extra strips should be immediately resealed with desiccant and returned to proper storage. Strips are stable for 60 days after the envelope has been opened and properly resealed and the indicator strip on the desiccant pouch remains blue.
3. Conjugate: Store between 2° and 8 °C. DO NOT FREEZE.
4. Calibrator, Positive Control and Negative Control: Store between 2 ° and 8 °C.
5. TMB: Store between 2 ° and 8 °C.
6. Wash Buffer concentrate (10X): Store between 2 ° and 25 °C. Dilute wash buffer (1X) is stable at room temperature (20 ° to 25 °C) for up to 7 days or for 30 days between 2 ° and 8 °C.
7. Sample Diluent. Store between 2 ° and 8 °C.
8. Stop Solution: Store between 2 ° and 8 °C.

REFERENCES

1. Rapp CE and Hewetson JF: Infectious mononucleosis and the Epstein-Barr virus. *Am. J. Dis. Child.* 132:78-86, 1978.
2. Biggar RJ, Henle W, Fleisher G, Bocker J, Lennette ET, and Henle G: Primary Epstein-Barr virus infections in African infants; I. Decline of maternal antibodies at time of infection. *Int. J. Cancer* 22:239-243, 1978.
3. Fry J: Infectious Mononucleosis: Some new observations from a 15 year study. *J. Family Pract.* 10:1087-1089, 1980.
4. Lennette ET: Epstein-Barr virus IN: *Manual of Clinical Microbiology*, 46th edition. Lennette EH, Balows A, Hausler WJ, Shadomy HJ (eds). American Society for Microbiology, Washington, DC, 1985.
5. Fleisher G, Henle W, Henle G, Lennette ET, and Biggar RJ: Primary infection with Epstein-Barr virus in infants in the United States: Clinical and Serologic observations. *J. Infect. Dis.* 139:553-558, 1979.
6. Merlin TL: Chronic mononucleosis: Pitfalls in the laboratory diagnosis. *Human Path.* 17:2-8, 1986.
7. Sixbey JW, Nedrud JG, Raab-Traub N, Hanes RA, and Pagano JS: Epstein-Barr virus replication in oropharyngeal epithelial cells. *New Engl. J. Med.* 310:1225-1230, 1984.
8. Chang RS, Lewis JP, Reynolds RD, Sullivan MJ, and Neuman J: Oropharyngeal excretion of Epstein-Barr virus by patients with lymphoproliferative disorders and by recipients of renal homografts. *Annals of Intern. Med.* 88:34-40, 1978.



9. Jones JF, Ray G, Minnich LL, Hicks MJ, Kibler R, and Lucas DO: Evidence of active Epstein-Barr virus infection in patients with persistent, unexplained illness: Elevated anti-early antigen antibodies. *Annals of Intern. Med.* 102:1-6, 1985.

10. Dubois RE, Seeley JK, Brus I, Sakamoto K, Ballow M, Harada S, Bechtold TA, Pearson G, and Purlito DT: Chronic mononucleosis syndrome. *South Med. J.* 77:1376-1382, 1984.

11. Straus SE, Tosato G, Armstrong G, Lawley T, Preble OT, Henle W, Davey R, Pearson G, Epstein J, Brus I, and Blaese M: Persisting illness and fatigue in adults with evidence of Epstein-Barr virus infection. *Annals of Intern. Med.* 102:7-16, 1985.

12. Tobi M, Ravid Z, Feldman-Weiss V, Ben-Chetrit E, Morag A, Chowers I, Michaeli V, Shalit M, and Knobler H: Prolonged atypical illness associated with serological evidence of persistent Epstein-Barr virus infection. *Lancet* 1:61-63, 1982.

13. Evans AS, Neiderman JC, Cenabre LC, West B, and Richards VA: A prospective evaluation of heterophile and Epstein-Barr virus-specific IgM antibody tests in clinical and subclinical infectious mononucleosis: Specificity and sensitivity of the tests and persistence of antibody. *J. Infect. Dis.* 132:546-554, 1975.

14. Henle W, Henle GE, and Horwitz CA: Epstein-Barr virus specific diagnostic tests in infectious mononucleosis. *Human Path.* 5:551-565, 1974.

15. Reedman BM, and Klein G: Cellular localization of an Epstein-Barr virus (EBV)-associated complement-fixing antigen in producer and non-producer lymphoblastoid cell lines. *Int. J. Cancer* 11:499-520, 1973.

16. Henle G, Henle W, and Horwitz CA: Antibodies to Epstein-Barr virus-associated nuclear antigen in infectious mononucleosis. *J. Infect. Dis.* 130:231-239, 1974.

17. Henle G, Henle W, and Klein G: Demonstration of two distinct components in the early antigen complex of Epstein-Barr virus-infected cells. *Int. J. Cancer* 8:272-282, 1971.

18. Henle W, Ho H-C, Henle G, and Kwan HC: Antibodies to Epstein-Barr virus-related antigens in nasopharyngeal carcinoma. Comparison of active cases with long-term survivors. *J. Nat'l. Cancer Inst.* 51:361-369, 1973.

19. Horwitz CA, Henle W, Henle G, and Schmitz H: Clinical evaluation of patients with infectious mononucleosis and development of antibodies to the R component of the Epstein-Barr virus-induced early antigen complex. *Am. J. Med.* 58:330-338, 1975.

20. Sumaya CV: Endogenous reactivation of Epstein-Barr virus infections. *J. Infect. Dis.* 135:374-379, 1977.

21. Joncas J, Lapointe N, Gervais F, Leyritz M, and Wills A: Unusual prevalence of antibodies to Epstein-Barr virus early antigen in ataxia telangiectasia. *Lancet* 1:1160, 1977.

22. Akaboshi I, Jamamoto J, Katsuki T, and Matsuda I: Unique pattern of Epstein-Barr virus specific antibodies in recurrent parotitis. *Lancet* 2:1049-1051, 1983.

23. Henle W, and Henle G: Epstein-Barr virus infection and pregnancy. *J. Infect. Dis.* 147:982-986, 1983.

24. Fleisher G, and Bolognese R: Persistent Epstein-Barr virus infection and pregnancy. *J. Infect. Dis.* 147:982-986, 1983.

25. Sumaya CV: Serological testing for Epstein-Barr virus; Developments in interpretation. *J. Infect. Dis.* 151:984-987, 1985.

26. Horwitz CA, Henle W, Henle G, Rudnick H, and Lutts E: Long-term serological follow-up of patients for Epstein-Barr virus after recovery from infectious mononucleosis. *J. Infect. Dis.* 151:1150-1153, 1985.

27. Chang RS, Thompson H, and Pomerantz S: Epstein-Barr virus infection in homosexual men with chronic persistent generalized lymphadenopathy. *J. Infect. Dis.* 151:459, 1985.

28. Lennette ET, and Henle W: Epstein-Barr virus infections: Clinical and Serologic features. *Lab Mgmt.* 23-28, June, 1987.

29. Procedures for the collection of diagnostic blood specimens by venipuncture - Second Edition; Approved Standard (1984).

Published by National Committee for Clinical Laboratory Standards.

30. Procedures for the Handling and Processing of Blood Specimens. NCCLS Document H18-A, Vol. 10, No. 12, Approved Guideline, 1990.

31. U.S. Department of Labor, Occupational Safety and Health Administration: Occupational Exposure to Bloodborne Pathogens, Final Rule. *Fed. Register* 56:64175-64182, 1991.

32. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guidelines - 4th Edition (2010). CLSI Document GP44-A4 (ISBN 1-56238-724-3). Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, PA 19087.

 ISO 13485 ISO 9001			
 Diagnostic Automation/Cortez Diagnostics, Inc. 21250 Califa Street, Suite 102 and 116, Woodland Hills, California 91367 USA			
Date Adopted	2016-12-5		
REF 1415-2	AccuDiag™ - EBV-EA IgG ELISA		
<table border="1"> <tr> <td>EC</td> <td>REP</td> </tr> </table>	EC	REP	CEpartner4U, Esdoornlaan 13, 3951DB Maarn. The Netherlands. www.cepartner4u.eu
EC	REP		
Revision Date: 2016-9-22			