) ( Di	agno	stic	A	utor	natio	N/	/Cc	rte	5 A	)íaç	ŢNO	sti	<b>(</b> 5,	Inc	•
	M N	U	Ν	0	D		Α	G	Ν	0	S	Т	(	C S	-
Accu Free ELIS REF 4	Diag™ e PSA 5A Kit 322-16 .abel 2°c	-8°C 🗡	796	Tests			prosta positiv deterr the pe those serum cancel would effect	itic hyp ve scree ninatior ercentag with b PSA in r screen subset s on car	erplasia ning re o of fre ge of f enign p conjun ing in s quently ncer de	e and p esults. <sup>5</sup> e PSA prostati ction w elected reduct tection	rostat A pote levels. A is lov ic hyp vith tot I men e unn rates. E	itis, lea ential s " <sup>3</sup> Preli wer in erplasi tal PSA with ela eccessa 6	ading to olution t minary s patients a. <sup>2</sup> Thus a. <sup>2</sup> Thus a. <sup>2</sup> can im evated t ary pros	a large p to this pr studies ha s with pr , the me prove sp otal seru tate biop	oercer oblen ave su ostate asure ecifici m PSA osies
Free P	SA ELIS	4					The f-	PSA ELI	SA tes	tisas	olid pl	hase t	wo-site i	immunoa	ssay.
Method	Imr	Enzym nunoso	ne Link orbent	ked t Assay			anoth used a	cional a er anti-F as the ti	ntibod 'SA mo racer, 1	y is coa noclona The f-PS	ated o al antil SA mol	on the body la lecules	surface beled wi	of the m ith horser t in the st	adish radish
Principle	S	andwic	h Con	nplex			serum	are "	sandw	iched"	betwo	een th	ne two	antibodi	es. F
Detection Range		0-10	ng/ml	L			forma	tion of	the co	ated a	ntibod	ly-antig	gen- anti	ibody-enz	yme
Sample		100 µ	L Seru	m			unbou	ind anti idasa ar	body-e	nzyme <sup>-</sup> bound	tracers	s are re	emoved l	by washir	ig. Th
Specificity		9	8.7%				reaction	on. The	e inter	nsity o	of the	color	forme	d is pro	oport
Sensitivity		0.05	ng/m	L			conce	ntratior	of f-PS	5A pres	ent in	the sar	nple.	F	
Incubation Time		160 1	ninute	es											
Shelf Life	12 m	anufac	turing	m the g date			SP	ECIN	IEN (	OLL	ΕΟΤΙ	ON A	AND F	<b>PREPA</b>	RAT
							1. E	Blood sh erum sh	iould b nould b	e drawr e separ	n using ated f	; standa rom th	ard venij e red blo	puncture ood cells a	techn as soc



### **INTENDED USE**

The DAI Free-PSA Enzyme Immunoassay test kit is intended for the quantitative determination of f-PSA in human serum.

# SIGNIFICANCE AND SUMMARY

Human Prostate Specific Antigen (PSA) is a 33 kD serine proteinase which, in human serum, is predominantly bound to alpha 1-antichymotrypsin (PSA-ACT) and alpha 2-macroglobulin (PSA-AMG). Trace amounts of alpha 1-antitrypsin and inter-alpha trypsin inhibitor bound to PSA can also be found. Any remaining PSA is in the free form (f-PSA).<sup>1-3</sup> Current methods of screening men for prostate cancer utilize the detection of the major PSA-ACT form. Levels of 4.0 ng/ml or higher are strong indicators of the possibility of prostatic cancer.<sup>4</sup> However, elevated serum PSA levels have also been attributed to benign

ge of false nvolves the gested that ancer than ent of free of prostate vels, which th minimal

anti-f-PSA wells and eroxidase is solution or owing the mplex, the orseradish olorimetric al to the

# **DN**

- ues and the as possible. Avoid grossly hemolytic, lipemic, or turbid samples.
- Plasma samples collected in tubes containing EDTA, heparin, or oxalate 2. may interfere with the test procedures and should be avoided.
- Specimens should be capped and may be stored up to 48 hours at 2-8°C, 3. prior to assaying. Specimens held for a longer time can be frozen at -20°C.
- Thawed samples must be mixed prior to testing. 4.

# MATERIALS AND COMPONENTS

#### Materials provided with the test kit

- Antibody-coated microtiter plate with 96 wells. 1.
- 2. Sample Diluent, 12 ml.
- з. Reference standards containing 0, 0.1, 0.5, 2.0, 5.0, and 10.0 ng/ml f-PSA, liquid (ready to use) or lyophilized form. 1 set.
- 4. Enzyme Conjugate Reagent, 22 ml.
- TMB Substrate, 12ml 5.
- 6. Stop Solution, 12 ml.
- Wash Buffer Concentrate (50X), 15 ml. 7.

#### Materials required but not provided

- Precision pipettes: 0.10, 0.20, and 1.0 ml. 1.
- Disposable pipette tips. 2.
- Distilled water. з.
- Vortex mixer or equivalent. 4.
- Absorbent paper or paper towels. 5.
- 6. Graph paper.
- A microtiter plate reader with a bandwidth of 10nm or less and an optical 7. density range of 0-2.5 OD or greater at 450nm.

Diagnostic Automation/Cortez Diagnostics, Inc.

21250 Califa St, Suite 102 and 116, Woodland Hills, CA 91367 USA Phone: 818-591-3030, Fax: 818-591-8383 Email: onestep@rapidtest.com Website: www.rapidtest.com

Diagnostic Automation/Cortez Diagnostics, Inc.

# I M U N O D I A G N O S T I C

#### **REAGENT PREPARATION**

- 1. All reagents should be brought to room temperature (18-22°C) and mixed by gently inverting or swirling prior to use. Do not induce foaming.
- If reference standards are lyophilized, reconstitute each standard with 0.5 ml distilled water. Allow the reconstituted material to stand for at least 20 minutes. Reconstituted standards should be sealed and stored at 2-8°C.
- 3. Dilute 1 volume of Wash Buffer (50x) with 49 volumes of distilled water. For example, Dilute 15 ml of Wash Buffer (50x) into distilled water to prepare 750 ml of washing buffer (1x). Mix well before use.

### ASSAY PROCEDURE

- 1. Secure the desired number of coated wells in the holder.
- 2. Dispense 100 $\mu l$  of standards, specimens, and controls into appropriate wells.
- 3. Dispense 100µl of sample diluent into each well.
- 4. Thoroughly mix for 10 seconds. It is very important to have a complete mixing in this step.
- 5. Incubate at 37°C for 60 minutes
- 6. Remove the incubation mixture by emptying plate contents into a suitable waste container.
- 7. Rinse and empty the microtiter wells 5 times with washing buffer (1X).
- 9. Strike the wells sharply onto absorbent paper or paper towels to remove all residual water droplets.
- 10. Dispense 200µl of Enzyme Conjugate Reagent into each well. Gently mix for 5 seconds.
- 11. Incubate at 37°C for 60 minutes.
- 12. Remove the incubation mixture by emptying plate contents into a suitable waste container.
- 13. Rinse and empty the microtiter wells 5 times with washing buffer (1X).
- 14. Strike the wells sharply onto absorbent paper to remove residual water droplets.
- 15. Dispense 100µl TMB solution into each well. Gently mix for 5 seconds.
- 16. Incubate at room temperature for 20 minutes in the dark.
- 17. Stop the reaction by adding 100 $\mu$ l of stop solution to each well.
- 18. Gently mix for 30 seconds to make sure that the blue color changes completely to yellow.
- 19. Using a microtiter plate reader, read the optical density at 450nm within 15 minutes.

#### **Important Notes:**

- 1. The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.
- 2. It is recommended that if manual pipetting is used, no more than 32 wells be used for each assay run, since pipetting of all standards, specimens and controls should be completed within 5 minutes. A full plate of 96 wells may be used if automated pipetting is available.
- 3. Duplication of all standards and specimens, although not required, is recommended.

### RESULTS

Calculate the mean absorbance value ( $A_{450}$ ) for each set of reference standards, controls, and patient samples. Construct a standard curve by plotting the mean absorbance obtained from each reference standard against its concentration in ng/ml on graph paper. The absorbance values are placed on the vertical, or Y-axis, and concentrations on the horizontal, or X-axis. Use the mean absorbance values for each specimen to determine the corresponding concentration of f-PSA in ng/ml from the standard curve.

### EXAMPLE OF STANDARD CURVE

Results of typical standard run with optical density reading at 450nm shown in the Y-axis against f-PSA concentrations shown in the X-axis.

This standard curve is for the purpose of illustration only, and should not be used to calculate unknowns. Each user should obtain his or her own standard curve and data.

f-PSA (ng/ml)	Absorbance (450nm)
0	0.002
0.1	0.065
0.5	0.262
2.0	0.962
5.0	1.976
10.0	2.991



### **EXPECTED VALUES AND SENSITIVITY**

From the literatures, clinical studies have shown that higher tota PSA levels and lower percentages of free PSA are associated with higher risks of prostate cancer. In one study, 52 individuals with benign prostate hyperplasia (BPH) and 77 individuals with prostate cancer were analyzed. The f-PSA/t- PSA ratio were summarized in the table below.

f-PSA/t-PSA Ratio						
	Median	Min	Max	Mean (95% confidence interval)		
BPH	0.18	0.04	0.42	0.19 (0.17-0.21)		
Prostate Cancer	0.09	0.02	0.53	0.12 (0.10-0.14)		

The relationship between f-PSA/t-PSA ratio and risk of prostate cancer is also age related. When total PSA is in the range of 4.0-10.0 ng/mL, a f-PSA/t-PSA ratio  $\leq$  0.10 indicates 49% to 65% risk of prostate cancer depending on age; a f-PSA/t-PSA ratio > 0.25 indicates a 9% to 16% risk of prostate cancer, depending on age.

Percent probability of prostate cancer							
f-PSA/t-PSA ratio	50-59 years	60-69 years	≥ 70 years				
≤ 0.10	49.2	57.5	64.5				
0.11-0.18	26.9	33.9	40.8				
0.19-0.25	18.3	23.9	29.7				
>0.25	9.1	12.2	15.8				

Diagnostic Automation/Cortez Diagnostics, Inc.

21250 Califa St, Suite 102 and 116, Woodland Hills, CA 91367 USA Phone: 818-591-3030, Fax: 818-591-8383 Email: <u>onestep@rapidtest.com</u> Website: <u>www.rapidtest.com</u>



Multiple factors such as population, age, specificity of test method may affect interpretation of f-PSA and t-PSA values. These ranges should be used as guidelines only. Each laboratory should establish its own reference values.

#### **Cross Reactivity**

Antigens	Concentration	% Cross-react.
PSA-ACT	500 ng/ml	0.2
AFP	10,000 ng/mL	0
CEA	5,000 ng/mL	0
CA 125	1,000 U/mL	0
CA 15-3	1,000 U/mL	0
CA 19-9	1,000 U/mL	0
∞-HCG	1,000 ng/ml	0
ß-HCG	1,000 ng/mL	0
HCG	50,000 mIU/ml	0

#### Precision

Intra-Assay			
	Replicates	<u>S.D.</u>	<u>% CV</u>
Level I	20	0.005	13.1
Level II	20	0.011	4.4
Level III	20	0.128	3.2

#### Linearity

Two patient sera were serially diluted with o ng/ml standard. The average recovery was 108.7%.

Sample A						
<b>Dilution</b>	Expected	Observed	<u>% Recovery</u>			
undiluted	8.691	8.691	100			
2X	4.346	4.455	102.5			
4X	2.173	2.290	105.4			
8X	1.086	1.258	115.8			
16X	0.543	0.617	113.6			
32X	0.272	0.294	108.1			
64X	0.136	0.147	108.1			
Average Recovery: 107.6%						

Sample B						
<b>Dilution</b>	Expected	<b>Observed</b>	<u>% Recovery</u>			
undiluted	7.015	7.015	100			
2X	3.508	3.516	100.2			
4X	1.754	1.834	104.6			
8X	0.877	0.970	110.6			
16X	0.438	0.509	116.2			
32X	0.219	0.249	113.7			
64X	0.110	0.136	123.6			
Average Recovery: 100.8%						

#### Recovery

Equal parts of diluted patient sera were mixed to test for interference by unknown materials, such as drugs or hormones, in the assay. Concentrations

of Free PSA were determined before (original and added) and after (observed). The average recovery was 99.2%.

	Sample 1						
<u>Sample</u>	<u>Orig.</u> <u>Conc</u>	<u>Added</u>	Expected	<u>Observed</u>	<u>% Recovery</u>		
А	9.802	0.141	4.972	4.540	91.3		
В	9.802	0.281	5.042	4.638	92.0		
С	4.699	2.168	3.434	3.342	97.3		
D	2.168	1.170	1.669	1.661	99.5		
E	0.563	0.281	0.422	0.423	100.2		
F	4.699	1.125	2.912	2.727	93.6		
G	0.563	0.141	0.352	0.399	113.4		
	Average Recovery: 98.2%						

Sample 2							
<u>Sample</u>	<u>Orig.</u> <u>Conc</u>	<u>Added</u>	Expected	<u>Observed</u>	<u>% Recovery</u>		
А	6.825	0.098	3.462	3.301	95.3		
В	6.825	0.195	3.510	3.156	89.9		
C	3.342	1.563	2.453	2.503	102.0		
D	1.706	0.781	1.244	1.205	96.9		
E	0.391	0.195	0.293	0.327	111.6		
F	3.342	0.781	2.062	1.964	95.2		
G	0.391	0.098	0.245	0.271	110.6		

Average Recovery: 100.2%

#### LIMITATIONS OF THE PROCEDURE

There are some limitations of the assay:

- 1. As with all diagnostic tests, a definite clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- 2. Studies have implicated possible interference in immunoassay results in some patients with known rheumatoid factor and antinuclear antibodies. Serum samples from patients who have received infusions containing mouse monoclonal antibodies for diagnostic or therapeutic purposes, may contain antibody to mouse protein (HAMA). Although we have added some agents to avoid the interferences, we cannot guarantee it will eliminate all the effects of that.

### STORAGE

- Unopened test kits should be stored at 2-8°C upon receipt and the microtiter plate should be kept in a sealed bag with desiccants to minimize exposure to damp air. The test kit may be used throughout the expiration date of the kit
- 2. (One year from the date of manufacture). Refer to the package label for the expiration date. Opened test kits will remain stable until the expiring date shown, provided it is stored as prescribed above.
- 2. A microtiter plate reader with a bandwidth of 10nm or less and an optical density range of 0-2 OD or greater at 450nm wavelength is acceptable for use in absorbance measurement.

Diagnostic Automation/Cortez Diagnostics, Inc.

21250 Califa St, Suite 102 and 116, Woodland Hills, CA 91367 USA Phone: 818-591-3030, Fax: 818-591-8383 Email: <u>onestep@rapidtest.com</u> Website: <u>www.rapidtest.com</u>



### REFERENCES

- Christensson A, Bjork T, Nilsson O, et al. Serum Prostate Specific Antigen Complexed to α-1-Antichymotrypsin As An Indicator of Prostate Cancer. J. of Urol. 150:100-105; 1993.
- Lilja H, Christensson A, Dahlen U, et al. Prostate-specific antigen in serum occurs predominantly in complex with α-1-antichymotrypsin. Clin Chem. 1991;37:1618-1625.
- Stenman U-H, Leinonen J, Alfthan H, Rannikko S, Tuhkanen K, Alfthan O. A complex between prostate-specific antigen and α-1-antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostate cancer:assay of the complex improves clinical sensitivity for cancer. Cancer Res. 1991;51:222-226.
- 4. Catalona WJ, Smith DS, Ratliff TL, Basler JW. Detection of organconfined prostate cancer is increased through prostate-specific antigenbased screening. JAMA. 1993;270:948-954.
- Stamey TA, Yang N Hay AR, McNeal JE, Freiha, FS, Redwine E. Prostatespecific antigen as a serum marker for adenocarcinoma of the prostate. N Engl J Med. 1987;317-:909-916.
- 6. Catalona et al. Percentage of Free Serum PSA and Prostate Cancer Detection. JAMA. 1995;274, No. 15: 1214-1220.
- Catalona WJ, Smith DS, Ratliff TL, et al. Measurement of prostatespecific antigen in serum as a screening test for prostate cancer. N Engl J Med. 1991;324:1156-1161. Erratum: N Engl J Med. 1991;325:1324.
- Carter HB, Pearson JD, Metter J, et al. Longitudinal evaluation of prostate-specific antigen levels in men with and without prostate disease. JAMA. 1992;267:2215-2220.
- Smith DS, Catalona WJ. Rate of change in serum prostate-specific antigen levels as a method for prostate cancer detection. J Urol. 1994;152:1163-1167.
- Benson MC, Whang IS, Pantuck A, et al. Prostate specific antigen density: a means of distinguishing benign prostatic hypertrophy and prostate cancer. J Urol. 1992; 147:815-816.
- 11. Catalona WJ, Hudson MA, Scardino PT, et al. Selection of optimal prostate specific antigen cut-offs for early detection of prostate cancer: receiver operating characteristic curves. J Urol. 1994;152:2037-2042.
- Smith DS, Catalona WJ. The nature of prostate cancer detected through prostate specific antigen-based screening. J Urol. 1994;152:1732-1736.

