



# Diagnostic Automation/Cortez Diagnostics, Inc.



## I M M U N O D I A G N O S T I C S

### AccuDiag™ Beta 2 Glycoprotein 1 IgA ELISA

REF 1494-11

IVD See External Label 2°C 96 Tests

Beta 2 Glycoprotein 1 IgA ELISA	
Principle	Indirect ELISA
Detection	Qualitative and Semi-Quantitative
Sample	5 µL serum/plasma
Incubation Time	~ 75 min.
Sensitivity	98%
Specificity	91%
Shelf Life	12 Months from the manufacturing date

#### PRODUCT FEATURES

- Very easy to use with little training
- Highly specific and consistent assay
- Provides accurate results quickly
- Reading of results both visually and as absorbance data

#### SUMMARY OF ASSAY PROCEDURE

Step	Room temperature (20-25°C)	Volume	Incubation time
1	Sample dilution 1:101 = 5 µl / 500 µl		
2	Washing buffer (3 times)	350 µl	
3	Diluted samples, controls & calibrators	100 µl	30 minutes
4	Washing buffer (3 times)	350 µl	
5	Enzyme conjugate	100 µl	30 minutes
6	Washing buffer (3 times)	350 µl	
7	TMB Chromogenic Substrate	100 µl	15 minutes

8	Stop solution	100 µl	
9	Reading OD 450 nm		

#### INTENDED USE

The DIAGNOSTIC AUTOMATION  $\beta_2$ GP1 IgA Enzyme-linked Immunosorbent Assay (ELISA) is intended for the detection and semiquantitative determination of IgA antibodies to  $\beta_2$ GP1 in human sera or plasma. The results of the assay are to be used as an aid in the diagnosis of certain autoimmune disease thrombotic disorders, anti-phospholipid syndrome, SLE or lupus-like disorders.

#### SIGNIFICANCE AND SUMMARY

Cardiolipin autoantibodies (ACA) are described for various autoimmune diseases. The presence of anti-cardiolipin antibodies in systemic lupus erythematosus (SLE) can be related to the development of thrombocytopenia, in gynecology they are supposed to cause intrauterine death or recurrent abortion. Furthermore, anti-cardiolipin antibodies have been found in some non-thrombotic neurological disorders like cerebrovascular insufficiency, cerebral ischemia or chorea and in myocardial infarction. (1)

Recent studies have shown that a 50kD serum cofactor is required for anticardiolipin antibodies, to bind to cardiolipin which has been coated onto plastic plates. The cofactor has been identified as  $\beta_2$ -glycoprotein 1 also termed apolipoprotein H.  $\beta_2$ GP1 has been known as an in vitro inhibitor of the intrinsic blood coagulation pathway, ADP-dependent aggregation, and prothrombinase activity of activated platelets. (2-7)

It has become apparent that anticardiolipin antibody from patients with anti-phospholipid syndrome (APS) recognize a modified  $\beta_2$ GP1 structure and not cardiolipin, native  $\beta_2$ GP1 or an epitope structurally defined by both cardiolipin and  $\beta_2$ GP1. (2-6)

Galli et al. (3) and Viard, et al. (8) reported that anti-cardiolipin antibody derived from SLE and APS were directed to the  $\beta_2$ GP1 molecule coated on polystyrene plates. Koike and Matsuura showed conclusively that  $\beta_2$ GP1 is indeed the antigen to which many anticardiolipin antibody patients are actually binding and furthermore showed that the phospholipid merely serves to link the  $\beta_2$ GP1 to the solid phase. (2-9)

$\beta_2$ GP1 autoantibodies are found in the immunoglobulin classes IgG, IgM and IgA. The determination of IgM antibodies is a valuable indicator in the diagnosis of beginning autoimmune disease, whereas IgG and/or IgA antibodies will be found in progressive stages of manifested autoimmune disorders. IgA antibodies are often associated with IgG antibodies. The determination of IgA antibodies seems to have a greater validity in thrombosis and fetal loss. (10). Indications for determination of anti  $\beta_2$ GP1 antibodies are: SLE, Thrombosis, Thrombocytopenia, Cerebral Ischemia, Chorea, Epilepsy, Recurrent Abortion and Intrauterine Death.

#### ASSAY PRINCIPLE

Purified  $\beta_2$ GP1 antigens are coated on the surface of microwells. Diluted patient serum or plasma, and calibrators, are added to the wells. The Anti  $\beta_2$ GP1 specific antibodies, if present, bind to the antigens. All unbound materials are washed away. After adding enzyme conjugate, it binds to the antibody-antigen complex. Excess enzyme conjugate is washed off, and TMB Chromogenic substrate is added. The enzyme conjugate catalytic reaction is stopped at a specific time. The intensity of the color generated is proportional to the amount of IgA specific antibodies in the sample. The

Diagnostic Automation/Cortez Diagnostics, Inc.

21250 Califa St, Suite 102 and 116, Woodland Hills, CA 91367 USA Phone: 818-591-3030, Fax: 818-591-8383

Email: [onestep@rapidtest.com](mailto:onestep@rapidtest.com) Website: [www.rapidtest.com](http://www.rapidtest.com)





# Diagnostic Automation/Cortez Diagnostics, Inc.



## IMMUNO DIAGNOSTICS

according to their own established procedures. The followings are a suggestive guideline.

- Negative: < 20 SAU
- Low positive: 20 – 40 SAU
- Moderate positive: 40 - 70 SAU
- High positive: > 70 SAU

A positive result suggests the possibility of certain autoimmune disease thrombotic disorders. A negative result indicates no β<sub>2</sub> GP1 IgA antibody or levels below the detection limit of the assay.

### PERFORMANCE CHARACTERISTICS

#### Sensitivity, specificity, and accuracy:

A total of 75 samples were assayed with the DIAGNOSTIC AUTOMATION ELISA β<sub>2</sub>GP1 IgA (X values) and with a reference ELISA (1) (Y values). The correlation equation is:

$$Y = 0.9221 X + 0.5522 \quad R^2 = 0.9375 \quad (n = 75)$$

DAI ELISA β <sub>2</sub> GP1 IgA	N	Reference ELISA (1)		Total
		N	P	
	32 (D)	1 (B)	33	
	3 (C)	39 (A)	42	
	35	40	75	

$$\text{Sensitivity} = A / (A+B) = 39 / (39 + 1) = 98 \%$$

$$\text{Specificity} = D / (C+D) = 32 / (3 + 32) = 91 \%$$

$$\text{Accuracy} = (A+D) / (A+B+C+D) = (39 + 32) / (39 + 1 + 3 + 32) = 71 / 75 = 95 \%$$

A second reference ELISA (2) was used to test 3 samples which reference ELISA (1) tested for negative and DIAGNOSTIC AUTOMATION ELISA tested for positive. The results are all positive for the 3 samples. The samples tested for positive with reference ELISA (1) and negative with DIAGNOSTIC AUTOMATION ELISA remains a positive result when it was tested with the reference ELISA (2).

### Precision:

Statistic for CV, mean and SD were calculated for each of three samples from the results of 8 determinations in a single run for intra-assay. Inter assay precision was calculated from the result of 8 determinations of 8 different runs.

Intra-assay	n	Mean SAU	SD	% CV
Serum A	8	14.9	0.35	2.38
Serum B	8	30.8	1.39	4.52
Serum C	8	58.9	0.99	1.68

Inter-assay	n	Mean SAU	SD	% CV
Serum A	8	15.6	0.38	2.4

- Galli, M., P. Comfurius, C. Maassen, H.C. Hemker, M.H. De Baets, P.J.C. Van Breda-Vriesman, T. Barbui, R.F.A. Zwaal, and E.M. Bevers. Anticardiolipin antibodies (ACA) directed not to cardiolipin but to a plasma protein cofactor. *Lancet*, 335: 1544 – 1547, 1990.

Serum B	8	31.2	1.42	4.55
Serum C	8	59.3	1.05	1.77

### INTERFERENCE AND CROSS-REACTIVITY

DIAGNOSTIC AUTOMATION β<sub>2</sub>GP1 IgA test does not cross-react with the following positive samples tested: Rubella, Toxo, CMV, H. pylori, Measles, Mumps, VZV, RF and HSV.

### LIMITATIONS OF THE ASSAY

- Diagnosis cannot be made on the basis of anti β<sub>2</sub> GP1 results alone. These results must be used in conjunction with information from clinical evaluation and other diagnostic procedure.
- The clinical significance of β<sub>2</sub> GP1 antibodies in diseases other than SLE is currently under investigation.
- When negative anti β<sub>2</sub> GP1 titers are found in the presence of clinical indications, a lupus anticoagulant, anti-cardiolipin or other additional testing is indicated.
- It is to be expected that some samples can be anti-cardiolipin positive yet anti β<sub>2</sub> GP1 negative. The anti β<sub>2</sub> GP1 test is a more specific marker of thrombotic risk. The anticardiolipin test can produce false positive results due to cross-reactivity with dsDNA or certain infectious disease antibodies.

### PRECAUTIONS

- Potential biohazardous materials:  
The calibrator and controls contain human source components, which have been tested and found nonreactive for Hepatitis B surface antigen as well as HIV antibody with FDA licensed reagents. However, as there is no test method that can offer complete assurance that HIV, Hepatitis B virus, or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2, as recommended in the Centers for Disease Control / National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories." 1984
- Do not pipette by mouth. Do not smoke, eat, or drink in the areas in which specimens or kit reagents are handled.
- The components in this kit are intended for use as an integral unit. The components of different lots should not be mixed.
- This product contains components preserved with sodium azide. Sodium azide may react with lead and copper plumbing to form explosive metal azide. On disposal, flush with a large volume of water.
- To prevent injury and chemical burns, avoid contact with skin and eyes or inhalation and ingestion of the following reagents: Enzyme conjugate, TMB chromogenic substrate and Stop solution.

### REFERENCES

- Hughes, G.R.V., Harris, E.N. and Gharavi, A.E. The Anticardiolipin Syndrome. *The Journal of Rheumatology* 1986; 13, 3: 486 – 489.
- McNeil, H.P., R.J. Simpson, C.N. Chesteman, and S.A. Krilis. Anti-phospholipid antibodies are directed against a complex antigen that includes a lipid-binding inhibitor of coagulation: β<sub>2</sub>-glycoprotein (apolipoprotein H), *Proc Natl Acad Sci USA*, 87: 4120 – 4124, 1990.
- Matsuura, E., Y. Igarashi, M. Fujimoto, K. Ichikawa, and T. Koike. Anticardiolipin cofactor (s) and differential diagnosis of autoimmune disease. *Lancet*, 336: 177 – 178, 1990.

Diagnostic Automation/Cortez Diagnostics, Inc.

21250 Califa St, Suite 102 and 116, Woodland Hills, CA 91367 USA Phone: 818-591-3030, Fax: 818-591-8383

Email: [onestep@rapidtest.com](mailto:onestep@rapidtest.com) Website: [www.rapidtest.com](http://www.rapidtest.com)



5. Koike, T. and E. Matsuura. What is the “true” antigen for anticardiolipin antibodies? *Lancet*, 337: 671 – 672, 1991.
6. Matsuura, E., Y. Igarashi, M. Fujimoto, K. Ichikawa, T. Suzuki, T. Sumida, T. Yasuda, and T. Koike. Heterogeneity of anticardiolipin antibodies defined by the anticardiolipin cofactor. *J Immunol*, 148: 3885 – 3891, 1992.
7. Nimph, J., H. Wurm, and G.M. Kostner.  $\beta_2$ -glycoprotein (apo-H) inhibits the release reaction of human platelets during ADP-induced aggregation. *Atherosclerosis*, 63: 109 – 114, 1987.
8. Viard, J.P., Z. Amoura, and J.F. Bach. Association of anti- $\beta_2$ -glycoprotein I antibodies with lupus-type circulating anticoagulant and thrombosis in systemic lupus erythematosus. *Am J Med*, 93: 181 – 186, 1992.
9. Tsutsumi A. et.al. Antibodies to  $\beta_2$ -glycoprotein1 and clinical manifestations in patients with Systemic Lupus Erythematosus. *Arthritis and Rheumatism*, 39:1466 –1474, 1996.
10. Lee, B., Silver, R.M. et. al. Immunoglobulin A anti-  $\beta_2$ GP1 antibodies in women who experience unexplained recurrent spontaneous abortion and unexplained fetal death. *Am J Obstet Gynecol*, 9: 748, 2001.

### MANUFACTURER AND BRAND DETAILS

**ISO 13485:2016**



ISO 13485  
Quality  
Management for  
Medical Devices  
CERTIFIED

 **Diagnostic Automation/Cortez Diagnostics, Inc.**  
21250 Califa Street, Suite 102 and 116,  
Woodland Hills, California 91367 USA

<b>Date Adopted</b>	2023-09
<b>Brand Name</b>	AccuDiag™
<b>REF</b> 1494-11	AccuDiag™ - Beta 2 Glycoprotein 1 IgA ELISA
<b>EC</b> <b>REP</b>	CEpartner4U, Esdoornlaan 13, 3951 DB Maarn, The Netherlands. <a href="http://www.cepartner4u.eu">www.cepartner4u.eu</a>

Revision Date: 2016-05-25