



### AccuDiag™ Beta 2 Glycoprotein 1 IgG ELISA

REF 1495-11

IVD See External Lab 2-8°C 96 Tests

Beta 2 Glycoprotein 1 IgG ELISA	
Principle	Indirect ELISA
Detection	Quantitative
Sample	5 µL serum/plasma
Incubation Time	~ 65 min.
Sensitivity	91%
Specificity	92%
Shelf Life	12 Months from the manufacturing date

#### PRODUCT FEATURES

- Very easy to use with little training
- Highly specific and consistent assay
- Provides accurate results quickly
- Reading of results both visually and as absorbance data

#### INTENDED USE

The DIAGNOSTIC AUTOMATION β<sub>2</sub>GP1 IgG Enzyme-linked Immunosorbent Assay (ELISA) is intended for the detection and semiquantitative determination of IgG antibodies to β<sub>2</sub>GP1 in human sera or plasma. The results of the assay are to be used as an aid in the diagnosis of certain autoimmune disease thrombotic disorders, anti-phospholipid syndrome, SLE or lupus-like disorders.

#### SIGNIFICANCE AND SUMMARY

Anti-Cardiolipin autoantibodies (ACA) are described for various autoimmune diseases. The presence of anti-cardiolipin antibodies in systemic lupus erythematosus (SLE) can be related to the development of thrombocytopenia, in gynaecology they are supposed to cause intrauterine death or recurrent abortion. Furthermore, anti-cardiolipin antibodies have been found in some non-thrombotic neurological disorders like cerebrovascular insufficiency, cerebral ischemia or chorea and in myocardial infarction. (1)

Recent studies have shown that a 50kD serum cofactor is required for anticardiolipin antibodies, to bind to cardiolipin which has been coated onto plastic plates. The cofactor has been identified as β<sub>2</sub>-glycoprotein 1 also termed apolipoprotein H. β<sub>2</sub>GP1 has been known as an in vitro inhibitor of the intrinsic blood coagulation pathway, ADP-dependent aggregation, and prothrombinase activity of activated platelets. (2-7) It has become apparent that anticardiolipin antibody from patients with anti-phospholipid syndrome (APS) recognize a modified β<sub>2</sub>GP1 structure and not cardiolipin, native β<sub>2</sub>GP1 or an epitope structurally defined by both cardiolipin and β<sub>2</sub>GP1. (2-6)

Galli et al. (3) and Viard, et al. (8) reported that anti-cardiolipin antibody derived from SLE and APS were directed to the β<sub>2</sub>GP1 molecule coated on polystyrene plates. Koike and Matsuura showed conclusively that β<sub>2</sub>GP1 is indeed the antigen to which many anticardiolipin antibody patients are actually binding and furthermore showed that the phospholipid merely serves to link the β<sub>2</sub>GP1 to the solid phase. (2-9)

Anti-β<sub>2</sub>GP1 autoantibodies are found in the immunoglobulin classes IgG, IgM and IgA. The determination of IgM antibodies is a valuable indicator in the diagnosis of beginning autoimmune disease, whereas IgG and/or IgA antibodies will be found in progressive stages of manifested autoimmune disorders. IgA antibodies are often associated with IgG antibodies. The determination of IgA antibodies seems to have a greater validity in thrombosis and fetal loss. (10) Indications for determination of anti β<sub>2</sub>GP1 antibodies are: SLE, Thrombosis, Thrombocytopenia, Cerebral Ischemia, Chorea, Epilepsy, Recurrent Abortion and Intrauterine Death.

#### ASSAY PRINCIPLE

Purified β<sub>2</sub>GP1 antigens are coated on the surface of microwells. Diluted patient serum or plasma, and calibrators, are added to the wells. The Anti β<sub>2</sub>GP1 specific antibodies, if present, bind to the antigens. All unbound materials are washed away. After adding enzyme conjugate, it binds to the antibody-antigen complex. Excess enzyme conjugate is washed off, and TMB Chromogenic substrate is added. The enzyme conjugate catalytic reaction is stopped at a specific time. The intensity of the color generated is proportional to the amount of IgG specific antibodies in the sample. The results are read by a microwell reader, and compared in a parallel manner with calibrators.

#### SPECIMEN COLLECTION & PREPARATION

1. Collect blood specimens and separate the serum.
2. Specimens may be refrigerated at 2 - 8°C for up to seven days or frozen for up to six months. Avoid repetitive freezing and thawing of serum sample.

#### MATERIALS AND COMPONENTS

##### Materials provided with the test kits

- |                                                                            |                |
|----------------------------------------------------------------------------|----------------|
| 1. Microwell strips: β <sub>2</sub> GP1 antigen coated wells.              | 12 x 8 wells   |
| 2. Sample Diluent: Yellow color solution                                   | 50 ml /bottle  |
| 3. Washing concentrate 50x.                                                | 15 ml / bottle |
| 4. TMB Chromogenic Substrate: Amber bottle.                                | 12 ml / bottle |
| 5. Enzyme conjugate:                                                       | 12 ml / bottle |
| 6. Calibrator set (1:101 prediluted): 6.3, 12.5, 25, 50, 100, and 200 SGU. |                |
| 1.0 ml /vial                                                               |                |
| 7. Control set (1:101 prediluted): Negative and Positive controls.         |                |
| Ranges are indicated on each label.                                        | 1.0 ml / vial  |
| 8. Stop solution                                                           | 12 ml / bottle |





## I M M U N O D I A G N O S T I C S

A total of 75 samples were assayed with the DIAGNOSTIC AUTOMATION ELISA  $\beta_2$ GP1 IgG (X values) and with a reference ELISA (1) (Y values). The correlation equation is

$$Y = 1.0327 X + 1.0057 \quad R^2 = 0.9815 \quad (n = 75)$$

Diagnostic automation ELISA $\beta_2$ GP1 IgG	Reference ELISA (1)			
		N	P	Total
	N	49 (D)	2 (B)	51
P	4 (C)	20 (A)	24	
<b>Total</b>	<b>53</b>	<b>22</b>	<b>75</b>	

Relative sensitivity =  $A / (A+B) = 20 / (20 + 2) = 91 \%$   
 Relative specificity =  $D / (C+D) = 49 / (4 + 49) = 92 \%$   
 Agreement =  $(A+D) / (A+B+C+D) = (20 + 49) / (20 + 2 + 4 + 49) = 69 / 75 = 92 \%$

Among 4 samples which reference ELISA (1) tested for negative and DIAGNOSTIC AUTOMATION ELISA tested for positive, all 4 samples gave positive results with a second reference ELISA (2) assay.

**Comparison with a reference Cardiolipin ELISA kit**  
 A total of 77 samples were assayed with the DIAGNOSTIC AUTOMATION ELISA  $\beta_2$ GP1 IgG (X values) and with a reference Cardiolipin (2) (Y values). The correlation equation is

$$Y = 1.092 X + 3.4339 \quad R^2 = 0.7954 \quad (n = 77)$$

Diagnostic automation ELISA $\beta_2$ GP1 IgG	Reference Cardiolipin (2)			
		N	P	Total
	N	24 (D)	27 (B)	51
P	2 (C)	24 (A)	26	
<b>Total</b>	<b>26</b>	<b>51</b>	<b>77</b>	

Relative sensitivity =  $A / (A+B) = 24 / (24 + 27) = 47 \%$   
 Relative specificity =  $D / (C+D) = 24 / (2 + 24) = 92 \%$   
 Agreement =  $(A+D) / (A+B+C+D) = (24 + 24) / 77 = 62 \%$

The relative sensitivity of the IgG  $\beta_2$  GP1 appear to be low in comparison to those of the IgG anti- cardiolipin. It is expected because the cardiolipin as an antigen it will recognize antibodies to  $\beta_2$  GP1, phospholipid as well as infectious antigens, (most of these were positive syphilis).

**Expected value**  
 145 serum specimens obtained from normal, asymptomatic blood donors were tested with the DIAGNOSTIC AUTOMATION- $\beta_2$  GP1 IgG test. The mean SGU =4, SD =3.

**Precision**  
 Statistic for CV, mean and SD were calculated for each of three samples from the results of 8 determinations in a single run for intra-assay. Inter assay precision was calculated from the result of 8 determinations of 8 different runs.

Intra-assay	n	Mean SGU	SD	% CV
Serum A	8	16.3	1.17	7.17
Serum B	8	33.8	1.25	3.68
Serum C	8	67.1	4.55	6.78

Inter-assay	n	Mean SGU	SD	% CV
Serum A	8	16.5	1.39	7.94
Serum B	8	35.9	2.17	6.04
Serum C	8	69.4	2.83	4.07

**INTERFERENCE AND CROSS-REACTIVITY**  
 DIAGNOSTIC AUTOMATION  $\beta_2$ GP1 IgA test does not cross-react with the following positive samples tested: Rubella, Toxo, CMV, H. pylori, Measles, Mumps, EBV VCA, H. pylori, RF and ANA.

### LIMITATIONS OF THE ASSAY

1. Diagnosis cannot be made on the basis of anti  $\beta_2$  GP1 results alone. These results must be used in conjunction with information from clinical evaluation and other diagnostic procedure.
2. The clinical significance of  $\beta_2$  GP1 antibodies in diseases other than SLE is currently under investigation.
3. When negative anti  $\beta_2$  GP1 titers are found in the presence of clinical indications, a lupus anticoagulant, anti-cardiolipin or other additional testing is indicated.
4. It is to be expected that some samples can be anti-cardiolipin positive yet anti  $\beta_2$  GP1 negative. The anti  $\beta_2$  GP1 test is a more specific marker of thrombotic risk. The anticardiolipin test can produce false positive results due to cross-reactivity with dsDNA or certain infectious disease antibodies.

### PRECAUTIONS

1. Potential biohazardous materials:  
The calibrator and controls contain human source components, which have been tested and found nonreactive for Hepatitis B surface antigen as well as HIV antibody with FDA licensed reagents. However, as there is no test method that can offer complete assurance that HIV, Hepatitis B virus, or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2, as recommended in the Centers for Disease Control / National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories." 1984
2. Do not pipette by mouth. Do not smoke, eat, or drink in the areas in which specimens or kit reagents are handled.
3. The components in this kit are intended for use as an integral unit. The Components of different lots should not be mixed.
4. This product contains components preserved with sodium azide. Sodium azide may react with lead and copper plumbing to form explosive metal azide. On disposal, flush with a large volume of water.
5. To prevent injury and chemical burns, avoid contact with skin and eyes or inhalation and ingestion of the following reagents: Enzyme conjugate, TMB chromogenic substrate and Stop solution.

### STORAGE CONDITIONS

1. Store the kit at 2 – 8°C.
2. Always keep microwells tightly sealed in pouch with desiccants. It is recommended to use up all wells within 4 weeks after initial opening of the pouch.
3. The reagents are stable until expiration of the kit.
4. Do not expose test reagents to heat, sun, or strong light during storage or usage.

### REFERENCES



- Hughes, G.R.V., Harris, E.N. and Gharavi, A.E. The Anticardiolipin Syndrome. *The Journal of Rheumatology* 1986; 13, 3: 486 – 489.
- McNeil, H.P., R.J. Simpson, C.N. Chesteman, and S.A. Krilis. Anti-phospholipid antibodies are directed against a complex antigen that includes a lipid-binding inhibitor of coagulation:  $\beta_2$ - glycoprotein (apolipoprotein H), *Proc Natl Acad Sci USA*, 87: 4120 – 4124, 1990.
- Galli, M., P. Comfurius, C. Maassen, H.C. Hemker, M.H. De Baets, P.J.C. Van Breda-Vriesman, T. Barbui, R.F.A. Zwaal, and E.M. Bevers. Anticardiolipin antibodies (ACA) directed not to cardiolipin but to a plasma protein cofactor. *Lancet*, 335: 1544 – 1547, 1990.
- Matsuura, E., Y. Igarashi, M. Fujimoto, K. Ichikawa, and T. Koike. Anticardiolipin cofactor (s) and differential diagnosis of autoimmune disease. *Lancet*, 336: 177 – 178, 1990.
- Koike, T. and E. Matsuura. What is the “true” antigen for anticardiolipin antibodies? *Lancet*, 337: 671 – 672, 1991.
- Matsuura, E., Y. Igarashi, M. Fujimoto, K. Ichikawa, T. Suzuki, T. Sumida, T. Yasuda, and T. Koike. Heterogeneity of anticardiolipin antibodies defined by the anticardiolipin cofactor. *J Immunol*, 148: 3885 – 3891, 1992.
- Nimph, J., H. Wurm, and G.M. Kostner.  $\beta_2$ -glycoprotein (apo-H) inhibits the release reaction of human platelets during ADP-induced aggregation. *Atherosclerosis*, 63: 109 – 114, 1987.
- Viard, J.P., Z. Amoura, and J.F. Bach. Association of anti- $\beta_2$ -glycoprotein I antibodies with lupus-type circulating anticoagulant and thrombosis in systemic lupus erythematosus. *Am J Med*, 93: 181 – 186, 1992.
- Tsutsumi A. et.al. Antibodies to  $\beta_2$ -glycoprotein<sub>1</sub> and clinical manifestations in patients with Systemic Lupus Erythematosus. *Arthritis and Rheumatism*, 39:1466 –1474, 1996.
- Lee, B., Siliver, R.M. et. al. Immunoglobulin A anti-  $\beta_2$ GP<sub>1</sub> antibodies in women who experience unexplained recurrent spontaneous abortion and unexplained fetal death. *Am J Obstet Gynecol*, 9: 748, 2001.

<b>ISO 13485:2016</b>			
			
ISO 13485 Quality Management for Medical Devices <b>CERTIFIED</b>			
 <b>Diagnostic Automation/Cortez Diagnostics, Inc.</b> 21250 Califa Street, Suite 102 and 116, Woodland Hills, California 91367 USA			
<b>Date Adopted</b>	<b>2023-09</b>		
<b>Brand Name</b>	<b>AccuDiag™</b>		
<b>REF</b> 1495-11	<b>AccuDiag™- Beta 2 Glycoprotein 1 IgG ELISA</b>		
<table border="1"> <tr> <td>EC</td> <td>REP</td> </tr> </table>	EC	REP	<b>CEpartner4U, Esdoornlaan 13, 3951DB          Maarn. The Netherlands.</b> <a href="http://www.cepartner4u.eu">www.cepartner4u.eu</a>
EC	REP		
<b>Revision Date: 2016-05-25</b>			

### SUMMARY OF ASSAY PROCEDURE

Step	(20-25°C Room temp.)	Volume	Incubation time
1	Sample dilution 1:10 = 5 $\mu$ l / 500 $\mu$ l		
2	Washing buffer (3 times)	350 $\mu$ l	
3	Diluted samples, controls & calibrators	100 $\mu$ l	30 minutes
4	Washing buffer (3 times)	350 $\mu$ l	
5	Enzyme conjugate	100 $\mu$ l	30 minutes
6	Washing buffer (3 times)	350 $\mu$ l	
7	TMB Chromogenic Substrate	100 $\mu$ l	15 minutes
8	Stop solution	100 $\mu$ l	
9	Reading OD 450 nm		

### MANUFACTURER AND BRAND DETAILS

Diagnostic Automation/Cortez Diagnostics, Inc.  
 21250 Califa St, Suite 102 and 116, Woodland Hills, CA 91367 USA Phone: 818-591-3030, Fax: 818-591-8383  
 Email: [onestep@rapidtest.com](mailto:onestep@rapidtest.com) Website: [www.rapidtest.com](http://www.rapidtest.com)