



CORTEZ DIAGNOSTICS INC.

23961 Craftsman Road, Suite D/E/F, Calabasas, CA 91302

Tel: (818) 591-3030 Fax: (818) 591-8383

onestep@rapidtest.com

technicalsupport@rapidtest.com

www.rapidtest.com



See external label



2°C-30°C



Σ=96 tests



Cat # Multiple - Drug

DOA PANEL TEST

Multiple Drug

RapiDip™ InstaTest

(10 Any Combination)

FOR THE QUALITATIVE ASSESSMENT OF DRUGS AND THEIR METABOLITES IN HUMAN URINE

For in vitro Diagnostic and Forensic Use

INTENDED USE

All of DOA Panel Test is an immunochromatography based one step in vitro test. It is designed for qualitative determination of drug substances in human urine specimens. This assay may be used in the point of care setting. Below is a list of cut-off concentrations for each drug using our test.

Amphetamine	1000 ng/ml of d-amphetamine
Barbiturate	300 ng/ml of secobarbital
Benzodiazepine	300 ng/ml of oxazepam
Buprenorphine	10 ng/ml of Buprenorphine-3- β-d-glucuronide
Cocaine	300 ng/ml of benzoylecgonine
EDDP	100 ng/ml of EDDP
MDMA (Ecstasy)	500 ng/ml of MDMA
Methadone	300 ng/ml of methadone
Methamphetamine	1000 ng/ml of (+)methamphetamine
Opiate*	300 ng/ml of morphine
Opiate II*	2000 ng/ml of morphine
Oxycodone	100 ng/ml of oxycodone
Phencyclidine	25 ng/ml of phencyclidine
Tricyclic antidepressants	1000 ng/ml of Nortriptyline
Cannabinoid (THC)	50 ng/ml of 11-nor-Δ ⁹ -THC-9-COOH

This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/ mass spectrometry (GC/MS) has been established as the preferred confirmatory method by the Substance Abuse Mental Health Services Administration. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

Cut-off concentration of 2000 ng/ml for Opiates Test

SUMMARY AND EXPLANATION

Amphetamines are a class of potent sympathomimetic agents with therapeutic applications. The most common amphetamines are d-amphetamine and d,l-amphetamine. Amphetamines are central nervous stimulants that cause the neurotransmitters epinephrine, norepinephrine and dopamine to be released into the brain and body giving users feelings of euphoria, alertness, and increased energy. Chronic abuse of amphetamine leads to tolerance and drug reinforcement effect. Cardiovascular responses to amphetamine include increased blood pressure and cardiac arrhythmias. More acute responses produce anxiety, paranoia, hallucinations and psychotic behavior. Amphetamine is metabolized by a number of pathways. In general, acid urine promotes excretion whereas alkaline urine retards it. In 24 hours, approximately 79% of the amphetamine dose is excreted in acid urine and about 45% in alkaline urine. Typically, about 20% is excreted as unchanged amphetamine. Unchanged amphetamine can be detected up to 1 –2 days after use.

Barbiturates are a group of prescription drugs that are frequently abused. They can depress the central nervous system. Acute higher dose induces exhilaration, sedation and respiratory depression. More acute responses produce respiratory collapse and coma. The effects of short-acting barbiturates, such as secobarbital last 3 to 6 hours. The effects of long-acting barbiturates such as phenobarbital last 10 to 20 hours. Short-acting barbiturates normally remain detectable in urine for 4 to 6 days, while long-acting barbiturates can be detected for up to 30 days. Barbiturates are excreted in the urine in unchanged forms, hydroxylated derivatives, carboxylated derivatives and glucuronide conjugates.

Benzodiazepines are a class of widely prescribed central nervous system depressants which have anxiolytic, hypnotic, anticonvulsant and muscle relaxant effects. Chronic abuse can result in addiction and tardive dyskinesia. Acute higher doses lead to drowsiness, dizziness, muscle relaxation, lethargy, coma and possible death. The effects of benzodiazepines use last 4 – 8 hours. Many of the benzodiazepines share a common metabolic route, and are excreted as oxazepam and its glucuronide in urine. Oxazepam is detectable in the urine for up to 7 days after drug use.

Buprenorphine a derivative of thebaine, is an opioid that has been marketed in the United States as the Schedule V parenteral analgesic Buprenex. In 2003, based on a reevaluation of available evidence regarding the potential for abuse, addiction, and side effect, DEA reclassified buprenorphine from a Schedule V to a Schedule III narcotic. Buprenorphine resembles morphine structurally but has a longer duration of action than morphine and can be administered sublingually as an analgesic. In October 2002, FDA approved the use of a buprenorphine monotherapy product, Subutex, and a buprenorphine/naloxone combination product, Suboxone, for the treatment of opioid addiction. Subutex and Suboxone are the first narcotic drugs available under the US Drug Act (DATA) of 2003 for the treatment of opiate dependence that can be prescribed in the US in a physician's work place. It has also been shown that buprenorphine has abuse potential and may itself cause dependency. In addition, a number of deaths have been recorded as a result of overdose with intravenously injected buprenorphine in conjunction with other psychotropic drugs such as benzodiazepines. Buprenorphine is metabolized primarily by n-dealkylation to form glucuronide-buprenorphine and glucuronide-norbuprenorphine.

Cocaine Derived from the leaves of cocoa plant, cocaine is a potent central nervous system stimulant as well as a local anesthetic. Some of the psychological effects induced by cocaine are: euphoria,

confidence and a sense of increased energy, accompanied by increased heart rate, dilation of the pupils, fever, tremors and sweating. Continued ingestion of cocaine could induce tolerances and physiological dependency which leads to its abuse. Cocaine is used by smoking, intravenous, intranasal or oral administration and excreted in the urine primarily as benzoylecgonine in a short period. Benzoylecgonine has a biological half-life of 5 – 8 hours, which is much longer than that of cocaine (0.5 – 1.5 hours), and can be generally detected for 12 – 72 hours after cocaine use or exposure.

EDDP 2-Ethylidine-1,5-dimethyl-3,3-diphenylpyrrolidine, is the primary metabolite of methadone. Methadone is a controlled substance and is used for detoxification and maintenance of opiate dependant patients. Patients on methadone maintenance may exhibit methadone (parent) levels that account for 5-50% of the dosage and 3-25% of EDDP in urinary excretion during the first 24 hours. The detection of EDDP is more beneficial than traditional methadone screening, in that EDDP exists only in urine from individuals that ingested methadone. The tampering of specimens by spiking the urine with methadone can be prevented. Secondly, renal clearance of EDDP is not affected by urinary pH, therefore the EDDP test provides a more accurate result of methadone ingestion than the methadone parent screen.

MDMA Methylenedioxymethamphetamine (Ecstasy) is a designer drug first synthesized in 1914 by a German drug company for the treatment of obesity. Those who take the drug frequently report adverse effects, such as increased muscle tension and sweating. MDMA is not clearly a stimulant, although it has, in common with amphetamine drugs, a capacity to increase blood pressure and heart rate. MDMA does produce some perceptual changes in the form of increased sensitivity to light, difficulty in focusing, and blurred vision in some users. Its mechanism of action is thought to be via release of the neurotransmitter serotonin. MDMA may also release dopamine, although the general opinion is that this is a secondary effect of the drug. The most pervasive effect of MDMA, occurring in almost all people who have taken a reasonable dose of the drug, is to produce a clenching of the jaws. The MDMA Ecstasy Test Strip yields a positive result when Methylenedioxymethamphetamine in urine exceeds 500ng/ml.

Methadone is a synthetic opioid, clinically available. It is used clinically for the treatment of severe pain and in maintenance programs for morphine and heroin addicts. Methadone acts on the central nervous and cardiovascular systems to produce respiratory and circulatory depression. Methadone also produces miosis and increases the tone of smooth muscle in the lower gastrointestinal tract while decreasing the amplitude of contractions. Acute higher doses induce analgesia, sedation, respiratory depression and coma. After methadone administration, the major urinary excretion products are methadone and its metabolites, EDDP and EMDP. Large individual variations in the urine excretion of methadone are output of methadone from 5-22%. Typically, following a 5 mg oral dose, methadone and EDDP account for 5% of the dose in the 24-hour urine. In those individuals on maintenance therapy, methadone may account for 5 to 50% of the dose in the 24-hour urine and EDDP may account for 3 to 25% of the dose.

Methamphetamine is the most popular synthetic derivative of the amphetamines. It is a potent sympathomimetic agent with therapeutic applications. Acute large doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, reduced appetite, and a sense of increased energy and power. More acute response produces anxiety, paranoia, psychotic behavior, and cardiac dysrhythmias. Methamphetamine is excreted in the urine as amphetamine and oxidized and deaminated derivatives. However, 10-40% of methamphetamine is excreted unchanged. Methamphetamine is generally detectable in the urine for 3 to 5 days after use.

Opiate Opioid analgesics comprised of a large group of substances that control pain by depressing the central nervous system. Acute high dose used by abusers or addicts can cause depressed coordination, disrupted decision, decreased respiration, hypothermia and coma. Morphine is excreted unmetabolized and is the marker metabolic product of opiates. Morphine and morphine glucuronide is detectable in urine for several days after opiates dose.

Oxycodone is known as Oxycontin, Roxicodone and is an ingredient of Percodan, Percocet, Roxicet and Tylox. Oxycodone is a semi-synthetic opiates derived from opium. Like other opiates, oxycodone is characterized by its analgesic properties, and the tendency for users to form a physical dependency and develop tolerance with extended use. Oxycodone is usually administered in combination with non-opiate analgesics such as acetaminophen and salicylates for the relief of moderate to severe pain. Oxycodone is a central nervous system depressant that may cause drowsiness, dizziness, lethargy, weakness and confusion. Toxicity in an overdose of oxycodone can lead to stupor, coma, muscle flaccidity, severe respiratory depression, hypotension, and stripiac arrest. Oxycodone is metabolized by N- and O-demethylation. One of the metabolites, oxymorphone, is a potent narcotic analgesic, while the other, noroxycodone, is relatively inactive. Between 33 to 61% of a single dose of oxycodone is excreted in a 24 hour urine collection and consists of 13-19% free oxycodone, 7-29% glucuronide conjugated oxycodone, 13-14% glucuronide conjugated oxymorphone and an unknown amount of noroxycodone. The detection time window of oxycodone is 1-3 days following use.

Phencyclidine commonly known as PCP, is a hallucinogen which interacts with dopamine, cholinergic and adrenergic systems. It has dose dependent stimulant, depressant, hallucinogenic and psychological effects. PCP is mostly administered by oral or intravenously. Even moderate amount of PCP, from 5 to 100 ng/ml, can result in psychotic, violent and self-destruction. At high doses, from 100 to 500 ng/ml, PCP can cause convulsions, hypertension, prolonged coma, absent peripheral sensation, and even death. PCP is metabolized via hydroxylation, oxidation, and conjugation with glucuronic acid in the liver. About 10% of the dose is excreted in urine as unchanged drug. PCP can be detected in the urine for 7 to 8 days after drug administration. For chronic users, PCP may persist in urine for 2 to 4 weeks. The length of time following drug use for which a positive result may occur is dependent upon several factors, including the frequency and amount of drug, metabolic rate, excretion rate, drug half-life, and the drug user's age, weight, activity, and diet.

TCA Tricyclic antidepressants (TCAs) are type of prescription drugs for the treatment of depressive disorders. Tricyclic Antidepressants consists of two main chemical classes. The tertiary amines boost serotonin levels and are usually prescribed for insomnia, irritability and overstimulation; these include amitriptyline, Imipramine and doxepin. The secondary amines which include nortriptyline, desipramine and Protriptyline, enhance norepinephrine levels and are prescribed for fatigue; withdrawal and inertness. TCA abuse can result in respiratory depression. Convulsions, blood pressure deviation, severe cardiac conditions, and coma. TCAs are taken orally or sometimes by injection. TCAs are excreted in the urine mostly in the form of metabolites for up to ten (10) days.

THC The agents of Marijuana that cause various biological effects in humans are called cannabinoid. Cannabinoid is a central nervous stimulant that alters mood and sensory perceptions, produces loss of coordination, impairs short term memory, and produces symptoms of anxiety, paranoia, depression, confusion, hallucination, and increased heart rate. Large doses of cannabinoid could cause the development of tolerances and physiological dependency and lead to abuse. A tolerance to the cardiac and psychotropic effects can occur and withdrawal syndrome produces restlessness, insomnia, anorexia and nausea. Δ^9 -THC is the primary active ingredient in cannabinoids. The main metabolite excreted in the urine is 11-nor- Δ^9 -THC-9-COOH, which are found within hours of exposure and remain detectable in the urine for 3-10 days after smoking.

PRINCIPLE

Each component strip of DOA Panel is based on the principle of specific immunochemical reaction between antibodies and antigen to analyze particular compound in human urine specimen. The assay relies on the competition for binding antibody. When drug is present in the urine specimen, it competes with drug conjugate for the limited amount of antibody-dye conjugate. When the amount of drug is equal or more than the cut-off, it will prevent the binding of drug conjugate to the antibody. Therefore, a

positive urine specimen will not show a colored band on the test line zone, indicating a positive result, while the presence of a colored band indicates a negative result.

A control line is present in the test window to work as procedural control. This colored band should always appear on the control line zone if the test device is stored in good condition and the test is performed appropriately.

MATERIAL PROVIDED

1. DOA Panel Test device. The amount of each coated antigen and/or antibody on the strip is less than 1.0 mg for antigen conjugate and is less than 1.0 mg for goat anti-mouse IgG antibody.
Test zone: contains drug bovine protein antigen conjugates
Control zone: contains Goat anti-mouse IgG antibody
Conjugate pad: contains mice monoclonal anti-drug antibody.
2. Instruction for use.

MATERIAL REQUIRED BUT NOT PROVIDED

1. Urine collection container.
2. Timer or clock.

STORAGE AND STABILITY

The test device should be stored at 2 to 30°C and will be effective until the expiration date stated on the package. The product is humidity-sensitive and should be used immediately after being open. Any improperly sealed product should be discarded.

PRECAUTIONS

1. For in vitro diagnostic and forensic use only.
2. Do not use the product beyond the expiration date.
3. Handle all specimens as potentially infectious.
4. Humidity sensitive product, do not open foil pouch until it is ready to be tested.
5. Use a new urine specimen cup for each sample to avoid cross contamination.

SPECIMEN COLLECTION AND PREPARATION

Fresh urine does not require any special handling or pretreatment. Specimen should be collected in a clean, dry, plastic or glass container. If the assay is not performed immediately, urine specimen may be refrigerated at 2-8 °C or frozen up to 7 days. Specimens should be brought to room temperature before testing. Urine specimens exhibiting a large amount of precipitate or turbidity should be centrifuged or allowed to settle before testing. Avoid contact with skin by wearing gloves and proper laboratory attire.

QUALITY CONTROL

Good Laboratory practice recommends the daily use of control materials to validate the reliability of device. Control materials should be assayed as clinical specimen and challenging to the assay cutoff concentration, e.g., 25% above and below cutoff concentration. If control values do not fall within establish range, assay results are invalid. Control materials which are not provided with this test kit are commercially available.

The Cortez Drugs of Abuse Test provides a built-in process control with a different antigen/antibody reaction at the control region (C). This control line should always appear regardless the presence of drug or metabolite. If the control line does not appear, the test device should be discarded and the obtained result is invalid. The presence of this control band in the control region serve as 1) verification that sufficient volume is added, 2) that proper flow is obtained.

PROCEDURE

1. Bring all materials and specimens to room temperature.
2. Remove the test card from sealed foil pouch.
3. Place the sample pad end into the urine specimen being careful to hold each pad in the urine without touching the plastic card.
4. Hold the card in the urine for 10 seconds, remove from the urine and replace the cap.
5. Read the results at 5 ~ 10 minutes after adding the sample.

Do not interpret the result after 10 minutes.

INTERPRETATION OF RESULTS

Negative:

Two colored bands form on any strip of the card. The appearance of two colored bands, one in test line zone and the other in control line zone, indicates negative result for that particular test(s). The negative result does not indicate the absence of drug in the specimen; it only indicates the level of tested drug in the specimen is less than cut-off level.

Positive:

One colored band form on any strip of the card. One colored band appears in control line zone. No colored band is found in test line zone. This is an indication the level of tested drug(s) in the specimen is above the cut-off level.

Invalid:

If there is no colored band in control line zone of any strip, the test result is invalid. Retest the sample with a new device.

Note: A borderline(+/-) in test line zone should be considered negative result.

LIMITATION OF PROCEDURE

The assay is designed for use with human urine only. A positive result with any of the tests indicates only the presence of a drug/metabolite and does not indicate or measure intoxication. There is a possibility that technical or procedural error as well other substances in certain foods and medicines may interfere with the test and cause false results. Please refer "SPECIFICITY" section for lists of substances that will produce either positive results, or that do not interfere with test performance. If a drug/metabolite is found present in the urine specimen, the assay does not indicate frequency of drug use or distinguish between drug of abuse and certain foods and medicines.

EXPECTED RESULTS

The DOA Panel Test is a qualitative assay. It identifies the drug(s) in human urine at its cut-off concentration or higher. The concentration of the drug(s) can not be determined by this assay. The test is intended to distinguish negative result from presumptive positive result. All positive results must be confirmed using an alternate method, preferably GC/MS.

PERFORMANCE CHARACTERISTICS

A. Accuracy

The accuracy of the DOA test panels were evaluated in each component strip and in comparison to GC/MS method at the following concentration: d-amphetamine 1000ng/ml (AMP), secobarbital 300 ng/ml (BAR), oxazepam, 300 ng/ml (BZO), buprenorphine-3-β-d-glucuronide 10ng/ml (BUP), benzoylecgonine 300ng/ml (COC), EDDP 100ng/ml, methadone 300 ng/ml (MTD), (+)methamphetamine 1000 ng/ml (MET), Methylenedioxymethamphetamine 500 ng/ml (MDMA), phencyclidine 25 ng/ml (PCP), morphine 300 ng/ml (OPI), morphine 2000 ng/ml (OPI II), oxycodone 100ng/ml (OXY), and 11-nor-Δ⁹-THC-9-COOH 50ng/ml (THC). Nortriptyline 1000ng/ml (TCA).

The results of each component strip are listed below:

1. Amphetamine The accuracy of the amphetamine test was evaluated in comparison to GC/MS method at a cut-off of 1000 ng/ml. One hundred (100) urine specimens with GC/MS confirmed d-amphetamine concentration were evaluated in this study. The results are summarized and presented below:

CortezAMP Test	(-)		(+))		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between -25% and c/o)	Near cutoff positive (between c/o and +25%)	GC/MS Positive (greater than +25% cut off)	
Positive	1	4	5	33	90.5
Negative	47	6	2	2	91.4
Total	48	10	7	35	

Positive % agreement: 90.5, Negative % agreement: 91.4

Nine (9) specimens were found discrepant between the new screening method and the GC/MS method. When compared those data, 67% (6 out of 9) of the discrepancy specimens were found between +25% to -25% of cutoff concentration (750-1250 ng/ml).

2. Barbiturate The accuracy of the barbiturate test was evaluated in comparison to GC/MS at a cut-off of 300 ng/ml of secobarbital. One hundred and nineteen (119) urine specimens with GC/MS confirmed barbiturate concentration were evaluated in this study. The results are summarized and presented below:

CortezBAR Test	(-)		(+))		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between -25% and c/o)	Near cutoff positive (between c/o and +25%)	GC/MS Positive (greater than +25% cut off)	
Positive	0	0	8	36	100
Negative	59	13	3	0	95.8
Total	59	13	11	36	

Positive % agreement: 100, Negative % agreement: 95.8.

Three (3) specimens were found discrepant between the CortezBAR and GC/MS method. When compared those data, 100% (3 out of 3) of the discrepancy specimens were found between cut-off and +25% cut-off concentration (300 – 375 ng/ml).

3. Benzodiazepine The accuracy of the benzodiazepine test was evaluated in comparison to GC/MS at a cut-off of 300 ng/ml of oxazepam. One hundred and four (104) urine specimens with GC/MS confirmed oxazepam concentration were evaluated in this study. The results are summarized and presented below:

CortezBZD Test	(-)		(+))		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between -25% and c/o)	Near cutoff positive (between c/o and +25%)	GC/MS Positive (greater than +25% cut off)	
Positive	6	3	6	44	98.0
Negative	39	5	1	0	83.0
Total	45	8	7	44	

Positive % agreement: 98.0, Negative % agreement: 83.0

Ten specimens were found discrepant between the CortezBZD and GC/MS method. When compared those data, 40% (4 out of 10) of the discrepancy specimens were found between -25% and +25% cut-off concentration (225 – 375 ng/ml).

4. Buprenorphine The accuracy of the Buprenorphine test was evaluated in comparison to GC/MS at a cut-off of 10 ng/ml of buprenorphine-3-β-d-glucuronide. One hundred and one (101) urine specimens with confirmed buprenorphine-3-β-d-glucuronide concentrations were evaluated in this study. Borderline readings were recorded as negative. The results are summarized and presented below:

CortezBUP Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (less than -25% cut off)	Near cutoff negative (between – 25% and c/o	Near cutoff positive (between c/o and +25%	GC/MS Positive (greater than +25% cut off)	
Positive	0	1	12	34	97.9
Negative	42	8	1	0	98
Total	42	9	13	34	

Positive % agreement: 97.9, Negative % agreement: 98.

Two (2) specimens were found discrepant between the CortezBUP and GC/MS method. When compared those data, 50% (1out of 2) of the discrepancy specimens were found between -25% cut-off and cut-off concentration (7.5 – 10 ng/ml).

5. Cocaine The accuracy of the cocaine test was evaluated in comparison to GC/MS at a cut-off of 300 ng/ml of benzoylecgonine. One hundred and two (102) urine specimens with GC/MS confirmed benzoylecgonine concentration were evaluated in this study. The results are summarized and presented below:

CortezCOC Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between – 25% and c/o	Near cutoff positive (between c/o and +25%	GC/MS Positive (greater than +25% cut off)	
Positive	0	2	3	45	97.9
Negative	46	5	0	1	96.2
Total	46	7	3	46	

Positive % agreement: 97.9, Negative % agreement: 96.2

Three (3) specimens were found discrepant between the CortezCOC and GC/MS method. When compared those data, 67% (2 out of 3) of the discrepancy specimens were found between –25% and +25% cut-off concentration (225 – 300 ng/ml).

6. EDDP The accuracy of the methadone metabolite (EDDP) test was evaluated in comparison to GC/MS method at a cut-off of 100 ng/mL EDDP. One hundred and sixty (160) specimens with EDDP concentration confirmed by GC/MS were evaluated.

Cortez EDDP Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between –25% and c/o	Near cutoff positive (between c/o and +25%	GC/MS Positive (greater than +25% cut off)	
Positive	0	2	8	70	97.5
Negative	70	8	2	0	97.5
Total	70	10	10	70	

Positive agreement = 97.5% Negative agreement = 97.5%

Four (4) specimens were found discrepant between the Cortez EDDP and GC/MS method. When compared those data, 100% (4 out of 4) of the discrepancy specimens were found between –25% and +25% cut-off concentration (75 – 125 ng/ml).

7. MDMA The accuracy of the methamphetamine (MDMA) test was evaluated in comparison to GC/MS at a cut-off of 500 ng/ml of MDMA. One hundred and eleven (111) urine specimens with GC/MS confirmed MDMA concentration were evaluated in this study. The results are summarized and presented below:

Cortez MDMA Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between -25% and c/o	Near cutoff positive (between c/o and +25%	GC/MS Positive (greater than +25% cut off)	
Positive	0	2	7	45	94.5
Negative	46	8	3	0	96.4
Total	46	10	10	45	

Positive % agreement: 94.2, Negative % agreement: 96.3

Five (5) specimens were found discrepant between the RapidMDMA and GC/MS method. When compared those data, 100% (5 out of 5) of the discrepancy specimens were found between –25% and cut-off concentration (375 – 625 ng/ml).

8. Methadone The accuracy of the CortezMTD test was evaluated in comparison to GC/MS at a cut-off of 300 ng/ml of methadone. One hundred and nineteen urine specimens with confirmed methadone concentrations were evaluated in this study. The results are summarized and presented below:

Cortez MTD Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than –25% cut off)	Near cutoff negative (between –25% and c/o	Near cutoff positive (between c/o and +25%	GC/MS Positive (greater than +25% cut off)	
Positive	0	0	9	49	96.7
Negative	50	9	2	0	100
Total	50	9	11	49	

Positive % agreement: 96.7, Negative % agreement: 100.

Two specimens were found discrepant between the CortezMTD and GC/MS method. When compared those data, 100% (2 out of 2) of the discrepancy specimens were found between cut-off and +25% cut-off concentration (300 – 375 ng/ml).

9. Methamphetamine The accuracy of the methamphetamine test was evaluated in comparison to GC/MS at a cut-off of 1000 ng/ml of (+)methamphetamine. Ninety nine (99) urine specimens with GC/MS confirmed (+)methamphetamine concentration were evaluated in this study. The results are summarized and presented below:

Cortez MET Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than –25% cut off)	Near cutoff negative (between –25% and c/o	Near cutoff positive (between c/o and +25%	GC/MS Positive (greater than +25% cut off)	
Positive	0	0	4	41	91.8
Negative	44	6	4	0	100
Total	44	6	8	41	

Positive % agreement: 91.8, Negative % agreement: 100

Four (4) specimens were found discrepant between the CortezMET and GC/MS method. When compared those data, 100% (4 out of 4) of the discrepancy specimens were found between +25% cut-off concentration (1000 – 1250 ng/ml).

10. **Opiate** The accuracy of the opiates test was evaluated in comparison to GC/MS at a cut-off of 300 ng/ml of morphine. One hundred and twenty three urine specimens with GC/MS confirmed morphine and codeine concentrations were evaluated in this study. The results are summarized and presented below:

Cortez OPI Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between -25% and c/o)	Near cutoff positive (between c/o and +25%)	GC/MS Positive (greater than +25% cut off)	
Positive	0	2	5	70	94.9
Negative	35	7	1	3	95.5
Total	35	9	6	73	

Positive % agreement: 94.9, Negative % agreement: 95.5

Six (6) specimens were found discrepant between the CortezOPI and GC/MS method. When compared those data, 50% (3 out of 6) of the discrepancy specimens were found between -25% and +25% cut-off concentration (225 – 375 ng/ml).

11. **Opiate II** The accuracy of Cortez Opiates II Test Strip was evaluated in comparison to GC/MS at a cut-off of 2,000ng/ml of morphine. One hundred and eight urine specimens with confirmed morphine and codeine concentrations were evaluated in this study. Among them, there are 72 real samples and 36 diluted from these samples. The results are summarized and presented below:

CortezOPI II	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between -25% and c/o)	Near cutoff positive (between c/o and +25%)	GC/MS Positive (greater than +25% cut off)	
Positive	0	0	10	40	94.3
Negative	45	10	3	0	100.0
Total	45	10	13	40	

Positive % agreement: 94.3, Negative % agreement: 100.0.

Three specimens were found discrepant between Cortez Opiates Test Strip and GC/MS method. When compared those data, 100% (3 out of 3) of the discrepancy specimens were found between +25% cut-off concentration (2000 – 2500ng/ml).

12. **Oxycodone** The accuracy of the oxycodone test was evaluated in comparison to GC/MS method at a cut-off of 100 ng/ml. One hundred and forty (140) urine specimens with GC/MS confirmed oxycodone concentration were evaluated in this study. The results are summarized and presented below:

CortezOXY Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between -25% and c/o)	Near cutoff positive (between c/o and +25%)	GC/MS Positive (greater than +25% cut off)	
Positive	0	0	2	52	93.1
Negative	77	5	3	1	100
Total	77	5	5	53	

Positive % agreement: 93.1, Negative % agreement: 100

Four specimens were found discrepant between CortezOXY and the GC/MS method. When compared those data, 75% (3 out of 4) of the discrepancy specimens were found between cut-off and +25% of cutoff concentration (100-125 ng/ml).

13. **Phencyclidine** The accuracy of the PCP test was evaluated in comparison to GC/MS at a cut-off of 25 ng/ml of phencyclidine. Ninety nine urine specimens with GC/MS confirmed phencyclidine concentration were evaluated in this study. The results are summarized and presented below:

Cortez PCP Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between -25% and c/o)	Near cutoff positive (between c/o and +25%)	GC/MS Positive (greater than +25% cut off)	
Positive	1	0	5	36	87.2
Negative	45	6	2	4	98.1
Total	46	6	7	40	

Positive % agreement: 87.2, Negative % agreement: 98.1

Seven specimens were found discrepant between the CortezPCP and GC/MS method. When compared those data, 28.6% (2 out of 7) of the discrepancy specimens were found between -25% and +25% cut-off concentration (25-31.3 ng/ml).

14. **TCA** The accuracy of the TCA test was evaluated in comparison to GC/MS at a cut-off of 1000 ng/ml of Nortriptyline. One hundred (100) urine specimens with GC/MS confirmed Nortriptyline concentration were evaluated in this study. The results are summarized and presented below:

Cortez TCA Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between -25% and c/o)	Near cutoff positive (between c/o and +25%)	GC/MS Positive (greater than +25% cut off)	
Positive	0	1	5	39	93.6
Negative	48	4	2	1	98.1
Total	48	5	7	40	

Positive % agreement: 97.8, Negative % agreement: 94.5

Four specimens were found discrepant between the CortezTCA and GC/MS method. When compared those data, 75% (3 out of 4) of the discrepancy specimens were found between -25% and +25% cut-off concentration (750 – 1250 ng/ml).

15. **THC** The accuracy of the THC test was evaluated in comparison to GC/MS at a cut-off of 50 ng/ml of 11-nor- Δ^9 -THC-9-COOH. One hundred (100) and fourteen urine specimens with GC/MS confirmed 11-nor- Δ^9 -THC-9-COOH concentration were evaluated in this study. The results are summarized and presented below:

Cortez THC Test	(-)		(+)		Percent agreement with GC/MS
	GC/MS Negative (Less than -25% cut off)	Near cutoff negative (between -25% and c/o)	Near cutoff positive (between c/o and +25%)	GC/MS Positive (greater than +25% cut off)	
Positive	5	3	7	30	97.4
Negative	50	5	1	0	88.7
Total	54	8	8	30	

Positive % agreement: 84.1, Negative % agreement: 98.2

Eight (8) specimens were found discrepant between the CortezTHC and GC/MS method. When compared those data, 50% (4 out of 8) of the discrepancy specimens were found between –25% and +25% cut-off concentration (37.5 – 62.5 ng/ml).

B. Sensitivity

The cut-off concentrations (sensitivity level) of DOA panel test are determined to be: AMP 1000 ng/ml, BAR, 300 ng/ml, BZO 300 ng/ml, BUP 10 ng/ml, COC 300 ng/ml, EDDP 100 ng/ml, MDMA 500 ng/ml, MTD 300 ng/ml, MET 1000 ng/ml, OPI 300 ng/ml, OPI II 2000 ng/ml, OXY 100 ng/ml, PCP 25 ng/ml, TCA 1000 ng/ml and THC 50 ng/ml.

C. Precision

The precision of DOA panel tests were determined by conducting the test with spiked controls and interpreted the results by three individuals to verify the random error of visual interpretation. The results of 50% above and 50% below cut-off specimens are 100% agreed by three observers:

Tested Drug	Concentration (ng/ml)	Number tested	Corrected result	% Corrected result
AMP	500	40	40	100
	1500	40	40	100
BAR	150	40	40	100
	450	40	40	100
BZO	150	40	40	100
	450	40	40	100
BUP	5	40	40	100
	15	40	40	100
COC	150	40	40	100
	450	40	40	100
EDDP	50	40	40	100
	150	40	40	100
MTD	150	40	40	100
	450	40	40	100
MDMA	250	40	40	100
	750	40	40	100
MET	500	40	40	100
	1500	40	40	100
OPI I	150	40	40	100
	450	40	40	100
OPI II	1000	40	40	100
	3000	40	40	100
OXY	50	40	40	100
	150	40	40	100
PCP	12.5	40	40	100
	37.5	40	40	100
TCA	500	40	40	100
	1500	40	40	100
THC	25	40	40	100
	75	40	40	100

D. Specificity

The specificity for DOA panel test were tested by adding various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

1. Interference testing

The DOA panel test performance at cut-off level is not affected when pH and Specific Gravity ranges of urine specimen are at 4.5 to 9.0 and 1.005 to 1.035.

The following substances were tested and confirmed did not interfere with DOA panel tests at the listed concentrations.

Glucose	2000 mg/dl
Human albumin	2000 mg/dl
Human hemoglobin	10 mg/dl
Urea	4000 mg/dl
Uric acid	10 mg/dl

2. Specificity

The following table lists compounds that are detected by DOA panel test which produced positive results when tested at levels equal or greater than the concentrations listed below:

<u>Tests</u>	<u>Compounds</u>	<u>Cut-off (ng/ml)</u>
Amphetamine	D-Amphetamine	1,000
	D/L-Amphetamine	2,000
	(±)3,4Methylenedioxyamphetamine	2,500
	I-Amphetamine	30,000
	(+)methamphetamine	> 100 µg/ml
	(±)3,4Methylenedioxymethamphetamine	> 100 µg/ml
Barbiturate	Alphenal	100
	Barbital	150
	Pentobarbital	150
	Phenobarbital	150
	Amobarbital	300
	Secobarbital	300
	Butalbital	5,000
Benzodiazepines	Nitrazepam	100
	Chloradiazepoxide HCl	300
	Clobazam	300
	Desmethyldiazepam	300
	Oxazepam	300
	Temazepam	300
	Alprazolam	1000
	Bromazepam	1000
	Diazepam	1000
	Flunitrazepam	1000
	Lorazepam	1000
	Clonazepam	2000
	Flurazepam	100
Buprenorphine	Buprenorphine-3-β-d-glucoronide	10
	Buprenorphine	200
Cocaine	Benzoylcgonine	300
	Cocaine	30,000
EDDP	EDDP	100
	EMDP	200,000
	Methadone	500,000
Methadone	Methadone	300
	Methadol	300
	Norpropoxyphene	> 100,000

MDMA	(+)3,4Methylenedioxyamphetamine (Ecstasy)	500
	(+)3,4-MDA	500
	(+)Methamphetamine	100,000
Methamphetamine (Ecstasy)	(+)Methamphetamine	1000
	(±)3,4Methylenedioxyamphetamine (Ecstasy)	1000
	d-Amphetamine	> 100, 000
	l-Amphetamine	> 100, 000
	(±)3,4Methylenedioxyamphetamine	> 100, 000
	Chloroquine	> 100, 000
	(-)Ephedrine	> 100, 000
	β-Phenylethylamine	> 100, 000
	Procaine	> 100, 000
	d-Pseudoephedrine	> 100, 000
Ranitidine	> 100, 000	
Opiate	Morphine	300
	Morphine-3-β-glucuronide	300
	Codeine	300
	Ethylmorphine	300
	Hydromorphone	300
	Nalorphine	750
	Heroin	1250
	Hydrocodone	1250
	Normorphine	2000
	Norcodeine	2500
	Naloxone	25,000
	Natrexone	100,000
	Oxycodone	> 100 µg/ml
Opiate II	Ethylmorphine	1,000
	Morphine	2,000
	Morphine-3-β-glucuronide	2,000
	Codeine	2,000
	6-Acetylmorphine	2,000
	Dihydrocodone	2,000
	Heroin	5,000
	Hydrocodone	7,500
	Hydromorphone	7,500
	Nalorphine	15,000
	Normorphine	20,000
	Norcodeine	100,000
	Naloxone	100,000
Oxycodone	100,000	
Oxycodone	Oxycodone	100
	Dihydrocodeine	20,000
	Codeine	100,000
	Hydromorphone	100,000
	Morphine	> 100 µg/ml
	Acetylmorphine	> 100 µg/ml
	Buprenorphine	> 100 µg/ml
	Ethylmorphine	> 100 µg/ml
Phencyclidine	PCP	25
	Tramadol	50,000
	N-Demethyl-cis-tramadol	100, 000
	O-demethyl-cis-tramadol	> 100, 000

	Norpropoxyphene	> 100, 000
TCA	Nortriptyline	1000
	Protriptyline	1000
	Imipramine	1000
	Desipramine	1000
	Amitriptyline	1000
	Doxepin	1000
	Nordoxepin	1000
	Promazine	500
	Trimipramine	2000
	Perphenazine	> 100, 000
	Chlorpromazine	> 100, 000
	Clomipramine	> 100, 000
	THC	11-nor-Δ^9-THC-9-COOH
11-nor-Δ^8-THC-9-COOH		37.5
11-hydroxy-Δ^9-THC		5000
Δ^8-Tetrahydrocannabinol		15000
Δ^9-Tetrahydrocannabinol		25000

The following compounds show no cross-reactivity at concentration up to 100 μ g/ml unless specified.

Acetaminophen	4-Acetamidophenol	Acetylsalicylic acid	Amikacin
Amitriptyline	Arterenol	Aspartame	Ascorbic acid
Atrophine	Caffeine	Camphor	Chloroquine
Chlopheniramine	Cortisone	Deoxyephedrine	Dextromethorphan
Digitoxin	Digoxin	Diphenhydramine	Ecgonine
Ecgonine methyl ester	Ephedrine	Epinephrine	Gentisic acid
Guaiacol glycer ester	Histamine	Hydrochlorothiazide	Homatrophine
Imipramine	Ibuprofen	soproterenol	Ketamine
Lidocaine	Meperidine	Methaqualon	Methylphenidate
Neomycin	Niacinamide	Perphenazine	Penicillin G
Phenylethylamine- α	Phenylpropanolamine	Promethazine	Pseudoephedrine
Quinine antidine	Salicylic acid	Tetracycline	Tetrahydrozoline
Theophylline	Thioridazine	Trifluoperazine	Tryptophan
Tyramine			

REFERENCES

1. Urine testing for drugs of abuse, NIDA Research Monograph 73 (1986)
2. Steven B. Karch, Drugs of abuse hand book, CRC Press, 1st. Ed. (1998)
3. Ray H. Liu and Bruce A. Goldberger, Handbook of workplace drug testing, AACCC Press, Washington DC (1995)

Date Adopted	Reference No.
2005-07-03	DA-DOA PANEL TEST-2008

CORTEZ DIAGNOSTICS, INC.

23961 Craftsman Road, Suite D/E/F, Calabasas, CA 91302

Tel: (818) 591-3030 Fax: (818) 591-8383

ISO 13485-2003



Revision Date: 10/10/08