

**OneStep**  
**FSH Urine**  
**RapiDip™ InstaTest**

**REF 198219-1-20**


 See external Label
  4-30°C
  1 Test

**Sensitivity**

**25 mIU/ml**

**INTENDED USE**


Cortez Diagnostics OneStep FSH Urine RapiDip™ InstaTest is a test kit for the determination of FSH (human Follicular Stimulating Hormone) concentration in urine specimens.

**TEST PRINCIPLE**

Test results are read visually without the need for any instrument.


**ASSAY PROCEDURE**

1. Make sure the pouch and sample are at room temperature. When you are ready to begin testing, open the sealed pouch by tearing along the notch. Remove the test kit from the pouch and use it as soon as possible.
2. Immerse the strip into the urine with the arrow end pointing towards the urine. Do not immerse past the MAX (marker line). Take the strip out after 3 seconds and lay the strip flat on a clean, dry, non-absorbent surface (such as the mouth of the urine container).
3. Wait 10-20 minutes and read result. It is important that the background is clear before the result is read. Do not read results after more than 30 minutes.



**RESULTS**

- **Negative:** Only one color band appears on the control region or the test band appears but is lighter than the control band. The FSH concentration of the sample is below the detection sensitivity of 25 mIU/mL.
- **Positive:** If two color bands are visible and the test band is equal to or darker than the control band, the FSH is present in the sample at or above the detection sensitivity of 25 mIU/mL.
- **Invalid:** No visible band at all or no colored band appears on the control (C) region. A total absence of color in both regions is an indication of procedure error and/or that test reagent deterioration had occurred. Repeat with a new test kit. If the problem persists, discontinue using the test kit immediately and contact your local distributor.



**PERFORMANCE CHARACTERISTICS**

**Cross Reactivity**

The cross reactivity of FSH test kits was evaluated with FSH homologous hormone. Homologous hormones, HCG, TSH, LH were added to urine samples containing FSH at concentration of 0, 25 or 100 mIU/mL. No cross reactivity was observed in the study (see Table 1).

**Non-specific interference**

One-Step FSH test was checked for possible interference from visibly hemolyzed, lipemic and icteric samples. Human

hemoglobin, bilirubin or albumin was spiked into samples with different concentration of FSH and tested using un-spiked sample as controls. No significant interference was observed in 20 sample testing results that were either positive or negative for FSH. The results are shown in Table 2.

**Table 1 – Cross-reactivity study of One-step FSH test kit**  
**Urine samples spiked with homologous hormones**



FSH conc. in urine sample (mIU/mL)	Unspiked urine samples	LH 1000mIU/mL	HCG 1000mIU/mL	TSH 1000mIU/mL
0	--	-	-	-
	-	-	-	-
	-	-	-	-
25	+	+	+	+
	+	+	+	+
	+	+	+	+
100	+	+	+	+
	+	+	+	+
	+	+	+	+

**Table 2 – Non-specificity study of One-step FSH test kits**

Sample No	Unspiked Sample	Urine samples spiked with (mg/mL)			
		Hemoglobin 10	1	Bilirubin 0.06	Albumin 100
1	-	-	-	-	-
2	-	-	-	-	-
3	-	-	-	-	-
4	-	-	-	-	-
5	-	-	-	-	-
6	-	-	-	-	-
7	-	-	-	-	-
8	-	-	-	-	-
9	-	-	-	-	-
10	-	-	-	-	-
11	+	+	+	+	+
12	+	+	+	+	+
13	+	+	+	+	+
14	+	+	+	+	+
15	+	+	+	+	+
16	+	+	+	+	+
17	+	+	+	+	+
18	+	+	+	+	+
19	+	+	+	+	+
20	+	+	+	+	+

## LIMITATIONS OF PROCEDURE

1. As with all diagnostic tests, all results must be considered with other clinical information available to the physician. A definite clinical diagnosis should only be made by the physician after all clinical and laboratory findings have been evaluated.
2. Single measurement of FSH level may give misleading result due to the pulsatile nature of FSH secretion. Samples drawn at different times of the day may vary. Elevated levels of FSH should be confirmed by more than one FSH measurement.
3. The test results should not be affected by pain relievers, antibiotics and other common drugs. After using the pill or patch that contains hCG or LH may affect the test result and should not be taken while using this test kit. In addition, the test will not work properly if you are pregnant or menopausal.
4. A positive FSH test only measures that FSH level is at or greater than 25 mIU/mL (above average) and the transition of menopause is suggested. Menopause by definition is having gone 12 months without a menstrual cycle. A positive FSH test does not give information on exactly when menopause will be reached.
5. If a negative result is given while some of the common symptoms associated with menopause are experienced, consult with the physician.
6. Alow FSH level may be observed if there is a presence of a tumor in hypothalamus.

<b>ISO 13485</b> <b>ISO 9001</b>   <b>Diagnostic Automation/  Cortez Diagnostics, Inc.</b> <b>21250 Califa St, Suite 102 and 116,  Woodland Hills, California 91367 USA</b>	
<b>Date Adopted</b>	<b>2016-01-31</b>
<b>REF</b> 198219-1-20	<b>CORTEZ- OneStep  FSH Urine  RapiDip™ InstaTest</b>
<b>EC</b> <b>REP</b>	<b>CEpartner4U, Esdoornlaan 13,  3951DB Maarn. The Netherlands.</b> <a href="http://www.cepartner4u.eu">www.cepartner4u.eu</a>
Revision Date: 2003-10-27	

## PRECAUTIONS

1. For in vitro diagnostic use only.
2. Do not use test kit beyond expiry date.
3. The test device should not be reused.
4. Keep out of the reach of children.
5. The test kit can be stored at temperatures between 4 to 30°C in the sealed pouch to the date of expiration.
6. The test kit should be kept away from direct sunlight, moisture and heat.